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A situation you have experienced time and time again:

Everything seems well-prepared for the removal of an unspecified dual rod system. But then, you realize there is no fitting screwdriver for the set screw or that the tulip-head pedicle screw does not match your screwdriver...

Thanks to our Spine Easy Screw Removal System, this situation will not create an unpredictable delay that could postpone your subsequent interventions - regardless of whether you are dealing with intact or broken pedicle screws or with older systems that still employed Schanz screws. With our Spine Easy Screw Removal System, you have a firm grip even on the most complicated explantations.

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
**COMPACT AND CLEARLY ARRANGED**

The Spine Easy Screw Removal System is your safe companion for the efficient removal of all dual rod instrumentation in the spine. This compact solution meets all the requirements of a modern and efficient explantation system. The intuitive design ensures that the system is self-explanatory and easy to use.

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HEX. Doppel-Bits, kanüliert  
 HEX. double-bit, cannulated

Größe / Dimension	Art.Nr. / Ref.No.
3.5 / 4.0 mm	92-ESR-HXHX3540
4.5 / 5.0 mm	92-ESR-HXHX3540
5.5 / 6.0 mm	92-ESR-HXHX3540

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**SMART UNIVERSAL PEDICLE SCREW REMOVAL COMPONENTS**

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**REMOVAL COMPONENTS FOR SCHANZ SCREWS, BROKEN AND COLD-WELDED SCREWS**

Größe / Dimension	Art.Nr. / Ref.No.
Ø 3.5 mm to 7.5 mm	92-UESR-KS3 to 92-UESR-KS6

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First select the appropriate bit for the set screw. To do this, connect the bit attachment that you think will fit to the quick coupling and place the bit on the screw head. Check that the selected size is right.

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Now, insert the bit into the screw head of the pedicle screw. Make sure that the slide ring is set to „counterclockwise“ for the ratchet function of the screw driver. Now carefully turn the screwdriver handle counterclockwise. Set screws are generally not particularly tight, and tend to be easy to loosen.

**REMOVING SET SCREWS**

All set screws can be easily removed with this system. Always make sure that you have selected the right HEX. or TX. bit. If none of the bits fit please use the KS extractor

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After selecting the correct bit, connect the quick coupling to the handle. You can choose between the T-handle and the I-handle.

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However, there are also exceptions, with very tight screws. In these situations, you must use the pedicle screw retrieval tool as a counter-torque wrench. Place the tool directly over the rod next to the screw you want to remove.

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To improve handling, you can also connect the bit quick coupling to the distance adapter. Now, turn the T-Handle counter-clockwise to loosen the set screw while continuing to apply counter-pressure on the rod. By applying counter-pressure, you avoid damaging the construction or tearing out adjacent screws.

## REMOVING SET SCREWS

All set screws can be easily removed with this system. Always make sure that you have selected the right HEX. or TX. bit. If none of the bits fit please use the KS extractor

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If none of the TX. or HEX. bits fit, you can also remove the screw plugs using the KS extractor. Please use the extractor that grips the sides of the screw head without touching its base. Connect the extractor directly to the T-handle.

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
If the internal screw heads are identical in size to the set screws, then use the same HEX. or TX. bit as you would to remove the set screws.

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If there is no fitting bit for the screw heads in the tulip, insert the rods with 4.0, 5.0 and 6.0 mm diameter. The rod inserts are placed in the tulip head with the rod pliers.

### REMOVING PEDICLE SCREWS

In many pedicle screw systems, the sizes of the screw heads that connect to the „tulip“ are identical to the sizes of set screws. However, pedicle screws can easily be removed with this clever system, even if they have a completely different fitting in the tulip-head.

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Connect the quick bit adapter to a handle and turn it counterclockwise to remove the pedicle screw.

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Now secure the rod insert with one of the previously removed set screws. Make sure that you tighten the set screw enough to firmly attach the pedicle screw to the tulip head.

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Now place the pedicle retrieval instrument over the tulip head and the rod.

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In some cases it may happen, that the set screw after removing the pedicle screw, still stuck and you can not remove the short rod. In this case, insert the pedicle screw with the tip ahead in the pedicle retrieval instrument. The rod fits into the recesses of this instrument.

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You can now easily remove the pedicle screw. Repeat this step until all screws are removed.

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Connect the pedicle retrieval instrument with the T-handle.


**SET-KONFIGURATION  
SPINE EASY SCREW  
REMOVAL SYSTEM**

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**SET-CONFIGURATION SPINE  
EASY SCREW REMOVAL  
SYSTEM**

Complete system available with the Ref. No.: 92-SET-S-ESR

Art.Nr. / Ref.No.	Beschreibung / Description	Menge Quantity
92-ESR-IN1	T-Handgriff m. Ratsche 1/4" Anschl. / T-handle w. ratchet, 1/4" conn.	(A) 1
92-ESR-IN3	Verlängerungsadapter mittel / Distance adapter medium	(B) 1
92-ESR-IN6	Bohrmaschinenadapter 1/4", für 3-Backen-Futter Adapter for drilling machine 1/4", triangular end	(C) 1
92-ESR-IN7	Schnellkupplung f. Doppel-Bits / Quick coupling f. double-bits	(D) 1
92-ESR-IN8	Extraktor für Spinal-Schrauben / Extractor f. spine screws	(E) 1
92-ESR-IN12	I-Griff m. Ratschenfunktion 1/4" Anschl. / I-Handle w. ratchet 1/4" conn.	(F) 1
92-ESR-IN13	Stabfasszange 45° gewinkelt / Rod grasping forceps 45° angled	(G) 1
92-ESR-IN9	Stabeinsatz 4.0 mm f. Extraktor / Rod inlay 4.0 mm f. extractor	(H) 2
92-ESR-IN10	Stabeinsatz 5.0 mm f. Extraktor / Rod inlay 5.0 mm f. extractor	(I) 2
92-ESR-IN11	Stabeinsatz 6.0 mm f. Extraktor / Rod inlay 6.0 mm f. extractor	(J) 2
92- ESR-HXHX3540	HEX.Doppel-Bit 3.5/4.0 mm / HEX.double-bit 3.5/4.0 mm	(K) 2
92- ESR-HXHX4550	HEX.Doppel-Bit 4.5/5.0 mm / HEX.double-bit 4.5/5.0 mm	(L) 2
92- ESR-HXHX5560	HEX.Doppel-Bit 5.5/6.0 mm / HEX.double-bit 5.5/6.0 mm	(M) 2
92-ESR-TXTX1520	TX.Doppel-Bit 15.0/20.0 mm / TX.double-bit 15.0/20.0 mm	(N) 2
92-ESR-TXTX2530	TX.Doppel-Bit 25.0/30.0 mm / TX.double-bit 25.0/30.0 mm	(O) 2
92-UESR-AS3	Extraktor f. 3.5-4.5 mm abgebr. Schrauben / Extractor f. 3.5-4.5 mm broken screws	(P) 2
92-UESR-AS4	Extraktor f. 4.5-5.5 mm abgebr. Schrauben / Extractor f. 4.5-5.5 mm broken screws	(Q) 2
92-UESR-AS5	Extraktor f. 5.5-6.5 mm abgebr. Schrauben / Extractor f. 5.5-6.5 mm broken screws	(R) 2
92-UESR-AS6	Extraktor f. 6.0-7.5 mm abgebr. Schrauben / Extractor f. 6.0-7.5 mm broken screws	(S) 2
92-UESR-KS3	Extraktor f. 3.5-4.5 mm kaltverschweißte Schrauben Extractor f. 3.5-4.5 mm cold welded screw	(T) 2
92-UESR-KS4	Extraktor f. 4.5-5.5 mm kaltverschweißte Schrauben Extractor f. 4.5-5.5 mm cold welded screw	(U) 2
92-UESR-KS5	Extraktor f. 5.5-6.5 mm kaltverschweißte Schrauben Extractor f. 5.5-6.5 mm cold welded screw	(V) 2
92-UESR-KS6	Extraktor f. 6.0-7.5 mm kaltverschweißte Schrauben Extractor f. 6.0-7.5 mm cold welded screw	(W) 2
92-TRAY-S-ESR-INO	Spinal ESR Instrumentensieb 480x250x60 mm Spine ESR instrument tray 480x250x60mm	1
92- CL-Standard1	Deckel f. Sieb 480x250 mm / Lid f. tray 480x250 mm	1

**Empfohlener Sterilisations-Container für 92-SET-S-ESR**  
**Recommended sterile container for 92-SET-S-ESR**

Art.Nr. / Ref.No.	Beschreibung / Description
99-110-100UESR	1/1 Sterilisations-Container, Deckel gelocht/Boden ungelocht, silber, außen 580x280x100 mm, innen 550x265x84 mm 1/1 sterile container, lid perforated/bottom non perforated, silver, outside 580x280x135 mm, inside 550x265x115 mm

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**THE EXTENSION FOR SPECIAL SYSTEMS**

Use our „extended“ system to remove add-on modules such as crossbars or older spinal implants that use „hex nuts.“ The extended system is the ideal supplement to the basic system.

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High end bit coupling with quick locking, secure grip, sleek design

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**THE RIGHT OPEN WRENCH FOR EVERY NUT**

Gabelschlüssel  
Double offset open wrench

Größe / Dimension
6 / 7 mm
8 / 9 mm

Gabelschlüssel  
Double offset open wrench

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
All open-end wrenches are angled at 45° for optimum use in situ.

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### EXTENDED SYSTEM

Complete system available with the  
Ref. No.: 92-SET-S-ESR-EXTENDED


Art.Nr. / Ref.No.	Beschreibung / Description		Menge / Quantity
92-ESR-IN7	Schnellkupplung f. Doppel-Bits / Quick coupling f. double-bits	(A)	1
92-ESR-HXHX1520	HEX.Doppel-Bit 1.5/2.0 mm / HEX.double-bit 1.5/2.0 mm	(B)	2
92-ESR-HXHX2530	HEX.Doppel-Bit 2.5/3.0 mm / HEX.double-bit 2.5/3.0 mm	(C)	2
92-ESR-TXTX0708	TX.Doppel-Bit 7.0/8.0 mm / TX.double-bit 7.0/8.0 mm	(D)	2
92-ESR-TXTX0910	TX.Doppel-Bit 9.0/10.0 mm / TX.double-bit 9.0/10.0 mm	(E)	2
92-ESR-TXHX2763	TX./HEX.Doppel-Bit 27/6.35 mm / TX./HEX.double-bit 27/6.35 mm	(F)	1
92-ESR-FLAT0810	Schlitz Doppel-Bit 0.8/1.0 mm / Slot double-bit 0.8/1.0mm	(G)	1
92-ESR-TSS6	Steckschlüssel 6 mm / Socket wrench 6 mm	(H)	2
92-ESR-TSS8	Steckschlüssel 8 mm / Socket wrench 8 mm	(I)	2
92-ESR-TSS9	Steckschlüssel 9 mm / Socket wrench 9 mm	(J)	2
92-ESR-TSS10	Steckschlüssel 10 mm / Socket wrench 10 mm	(K)	2
92-ESR-TSS11	Steckschlüssel 11 mm / Socket wrench 11 mm	(L)	2
92-ESR-TSS12	Steckschlüssel 12 mm / Socket wrench 12 mm	(M)	2
92-ESR-SSS0607	Doppelgabelschlüssel 6/7 mm, 45° abgewinkelt Double offset open wrench 6/7 mm, 45° angled	(N)	1
92-ESR-SSS0809	Doppelgabelschlüssel 8/9 mm, 45° abgewinkelt Double offset open wrench 8/9 mm, 45° angled	(O)	1
92-ESR-SSS1011	Doppelgabelschlüssel 10/11 mm, 45° abgewinkelt Double offset open wrench 10/11 mm, 45° angled	(P)	1
92-ESR-SSS1213	Doppelgabelschlüssel 12/13 mm, 45° abgewinkelt Double offset open wrench 12/13 mm, 45° angled	(Q)	1
92-TRAY-S-ESR-IN2	Spinal ESR Instrumentensieb, erweitert, 250x240x60 mm Spine ESR instrument tray, extended, 250x240x60 mm		1
92-CL-STANDARD1	Deckel f. Sieb 250x240 mm Lid f. tray 250x240 mm		1


### Empfohlener Sterilisations-Container für 92-SET-S-ESR-EXTENDED Recommended sterile container for 92-SET-S-ESR-EXTENDED

Art.Nr. / Ref.No.	Beschreibung / Description
99-310-135SDESR	1/2 Sterilisations-Container, Deckel gelocht/Boden ungelocht, silber, außen 285x285x135 mm, innen 262x260x116 mm 1/2 sterile container, lid perforated/bottom non perforated, silver, outside 285x280x135 mm, inside 262x260x116 mm

**DIGI**med   
Medizintechnik

Digi-med

Inhaber: Patrick Mohr

Kreutzerstrasse 1

78573 Wurmlingen

Telefon: +49-(0)7461-9101172

Fax: +49-(0)7461-9101172

info@digimed-de

www-digimed-de

Rev.: 2024