**DISTAL RADIUS
POLYAXIAL ANGLE STABLE**

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Digimed Medizintechnik – a committed specialist

that cares for clinical resources
Specialized in the fields of orthopedics, traumatology and spinal column, Digimed Medizintechnik develops innovative products and solutions, from which you as a dealer benefit just like the respectively responsible doctor or the OP personnel and - beyond the economic and administrative aspect - clinic management as well. The excellent quality and precision of our products, as well as the willingness to adapt our developments exclusively to your requirements, are the most important elements of our offered products and services. Our product ideas help you to preserve the most precious resources of daily clinic operations: Nerves, time and money.

Polyaxial Winkelstabil Implantate Kits für Frakturversorgungen	4
Polyaxial angle stable Implant kits for fracture treatments	4
Polyaxial Winkelstabil Schrauben	5
Polyaxial angle stable screws	5
Polyaxial Winkelstabil distale Radius Platten	6-7
Polyaxial angle stable distal radius plates	6-7
Polyaxial Winkelstabil single use Instrumenten-Kit	8
Polyaxial angle stable single use instrument kit	8
Wichtige Informationen des distalen Radius-Systems	9
Important information about the distal radius system	9
Polyaxiale Implantation der Polyaxial Winkelstabil distalen Radius Platte	10-13
Polyaxial implantation of the Polyaxial Winkelstabil distal radius plate	10-13
Monoaxiale Implantation der Polyaxial Winkelstabil distalen Radius Platte mit Eindrehborrhülse	14-17
Monoaxial implantation of the Polyaxial Winkelstabil distal radius plate width monoaxial drill sleeve	14-17
Retrax, Ihre geniale Hilfe bei allen Radius-Frakturen	18-20
Retrax, your clever support for all radius fractures	18-20
Polyaxial Winkelstabil distale Fibula Platten	21-22
Polyaxial angle stable distal fibula plates	21-22
Wichtige Informationen des distalen Fibula-Systems	23
Important information about the distal fibula system	23
Monoaxiale Implantation der Polyaxial Winkelstabil Fibula Platte	24-26
Monoaxial implantation of the Polyaxial Winkelstabil fibula plate	24-26
Polyaxial Winkelstabil single use Radius System - Empfohlene Grundausstattungsartikel	27
Polyaxial angle stable single use radius system - basic recommended items	27
Polyaxial Winkelstabil single use Fibula System - Empfohlene Grundausstattungsartikel	28
Polyaxial angle stable single use fibula system - basic recommended items	28
SET-RADIUS10 Polyaxial Winkelstabil Radius Instrumenten und Implantate Set, unsteril	29-30
SET-RADIUS10 Polyaxial angle stable radius instrument and implants set, non-sterile	29-30
SET-FIBULA10 Polyaxial Winkelstabil Fibula Instrumenten und Implantate Set, unsteril	31-32
SET-FIBULA10 Polyaxial angle stable fibula instrument and implants set, non-sterile	31-32

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Polyaxial angle stable IMPLANT KITS FOR FRACTURE TREATMENTS

The efficient and clever way to operate
Do you want to operate quickly and safely?
If yes, then Polyaxial Winkelstabil is the right product for you.

Safe because

» everything is always delivered correctly and in sterile packaging

» all processes from production to the sterilization of products are validated

» gap-free batch traceability is ensured at all times

» only the best production and implant materials are used (Material: Ti6Al-4V / titanium-anodized surface treatment, type II)

Quick because

» in-house preparations are no longer required

» everything delivered is already sterile

» all products are easy to implant

» after the operation, you are ready at once for the next procedure

Do you want to lower your own costs and use your staff more efficiently?

This is not a problem with Polyaxial angular stable, because

» we settle the costs of all implants directly with the insurance company or billing clearinghouse (within Germany)

» you no longer have to pay for the implants yourself (within Germany)

» we deliver and pick up all Polyaxial angle stable kits free of charge (within Germany)

» your staff does not have to spend long hours preparing instruments and implants

» your staff is relieved of many tedious administrative tasks

» in-house preparations are no longer required

» everything delivered is already sterile

» all products are easy to implant

» after the operation, you are ready at once for the next procedure

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GREAT BECAUSE IT'S SO EASY

Congratulations - you are using the innovative Polyaxial angle stable single use implant system.

Please note:

Each implant comes in its own sterile packaging and features stickers for patient documentation. Open and remove from the kit only the implants intended for the operation. Unpacked and contaminated implants may not be returned. The instruments must be disposed of properly after surgery. If necessary, we can provide you with special disposal containers.

All kits are intended only for a single operation. You only need the enclosed single use instruments kit when polyaxial screws are used. Use the separately packaged monoaxial single use screw-in drill sleeve for treatments with monoaxial screws.

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Dimensions:

- » Thread diameter
- » Core diameter
- » Screw head diameter
- » Screwdriver head
- » Different lengths in 2 mm increments

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* Color code only on sterile packaging.

Article also available non-sterile (Ref. No. 93-6127-010PTXL to 93-6127-036PTXL and Ref. No. 93-6127-010TX to 93-6127-050TX).

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Material:

- » Titanium alloy Ti-6Al-4V
 - » Surface treatment: type II titanium anodization
- Reduces cold welding*

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
Polyaxial angle stable DISTAL RADIUS PLATE


Perfect design and 100% Made in Germany.

Already 10 years ago, we redefined the shape of the radius plate. In contrast to conventional radius plates, the distal rim of the plate does not rise to the radial, but to the ulnar aspect. This special design allows the plate to be positioned very close to the joint space, thus achieving excellent support of the fragments. The Polyaxial angle stable radius plate imitates the anatomical bone surface. The fact that the distal screw holes are arranged on two rows ensures secure fixation of the fracture, even when dealing with multifragmentary fractures. The anatomical shape enhances stability and safety.


**POLYAXIAL WINKELSTABIL RADIUS PLATTEN LINKS
POLYAXIAL ANGLE STABLE RADIUS PLATE LEFT**


93-6600-203
Farbcode* / Color code*
Distale Radius Platte / distal radius plate,
7/3 Loch, kurz, links / 7/3 holes, short, left

93-6600-205
Farbcode* / Color code*
Distale Radius Platte
distal radius plate,
7/3 Loch, standard, links
7/3 holes, standard, left

93-6600-215
Farbcode* / Color code*
Distale Radius Platte
distal radius plate,
7/3 Loch, standard, breit, links
7/3 holes, standard, wide, left

93-6600-207
Farbcode* / Color code*
Distale Radius Platte / distal radius plate,
7/5 Loch, lang, links / 7/5 holes, long, left

* Farbcode nur auf Sterilverpackung.
Artikel auch unsteril erhältlich
(Art. Nr. 93-6600-203, 93-6600-205, 93-6600-215, 93-6600-207).
* Color code only on sterile packaging.
Article also available non-sterile
(Ref. No. 93-6600-203, 93-6600-205, 93-6600-215, 93-6600-207).

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

Polyaxial angle stable DISTAL RADIUS PLATE
Perfect design and 100% Made in Germany.


POLYAXIAL WINKELSTABIL RADIUS PLATTEN RECHTS
POLYAXIAL ANGLE STABLE RADIUS PLATE RIGHT


93-6600-202
Farbcode* / Color code*
Distale Radius Platte
distal radius plate
7/3 Loch, kurz, rechts
7/3 holes, short, right

93-6600-204
Farbcode* / Color code*
Distale Radius Platte
distal radius plate
7/3 Loch, standard, rechts
7/3 holes, standard, right

93-6600-214
Farbcode* / Color code*
Distale Radius Platte
distal radius plate
7/3 Loch, standard, breit, rechts
7/3 holes, standard, wide, right

93-6600-206
Farbcode* / Color code*
Distale Radius Platte
distal radius plate
7/5 Loch, lang, rechts
7/5 holes, long, right

* Farbcode nur auf Sterilverpackung.
Artikel auch unsteril erhältlich
(Art. Nr. 93-6600-202, 93-6600-204, 93-6600-214, 93-6600-206).
* Color code only on sterile packaging.
Article also available non-sterile
(Ref. No. 93-6600-202, 93-6600-204, 93-6600-214, 93-6600-206).

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Material:

- » Titanium alloy TI-6Al-4V
 - » Surface treatment: type II titanium anodization
- Reduces cold welding*

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**Polyaxial angle stable, EVERYTHING STERILE AND
SAFELY PACKED**

Polyaxial angle stable single use instrument kit

Polyaxial angle stable instrument kit:
A separate instrument kit is included in each Polyaxial
angular stable kit. So you can implant the Polyaxial
angle stable plates fast, simply and safely. The standard
instrument kit includes drill bit, drill sleeves, T10
screwdriver and Kirschner wire.

**POLYAXIAL WINKELSTABIL SINGLE USE INSTRUMENTEN-KIT
POLYAXIAL ANGLE STABLE SINGLE USE INSTRUMENT KIT**

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Kirschner wire with trocar point, round end,
1.6 mm/150 mm (2 pcs.)

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IMPORTANT INFORMATION YOU MUST KNOW

Indication

The 2.7 mm distal radius system is indicated for the treatment of intra- and extra- articular fractures and osteotomies of the distal radius.

Polyaxial stabilization

Using the polyaxial, funnel-shaped screw-in drill sleeve, you can insert the PTXL screws into the plate holes in a freely selectable angle of up to 15°. With polyaxial stabilization, it is advisable to check the direction of the drill with each drilling operation using an image converter to ensure that the joint surface is not damaged.

Monoaxial stabilization

Using the monoaxial screw-in drill sleeve, the PTXL screws are inserted into the plate holes in the direction determined by the plate design. With monoaxial stabilization, always fill the distal, second plate hole on the radial side first. A first check of the drilling process using the image converter is required only with this plate hole. If the drill channel is outside the joint surface, all other screws will not affect the joint surface either, which means monoaxial stabilization eliminates the constant monitoring with the image converter. Polyaxial Winkelstabil more than just plates and screws

Self-locking screws

The Polyaxial Winkelstabil T10 screwdriver is manufactured with very high precision. This high precision allows for the screwdriver tip being „wedged“ in the screw head so that it is firmly connected with the screwdriver. When positioning the screwdriver in the head of the screw, make sure to press the screwdriver down into the screw head using a certain amount of pressure. If using insufficient pressure, the tip of the screwdriver may not penetrate deep enough into the screw head. If this is the case, the screw head may become damaged during tightening.

The Polyaxial Winkelstabil implants are all packaged individually and sterile and are color-coded. Please note the correct side when selecting the plate. Check packaging and labeling on the plate to ensure that the plate is for the correct side („R“ for right and „L“ for left). The distal lip on the ulnar side of the radius plate is always the prominent one.

Soft tissue procedure

Make an approx. 5 cm long longitudinal incision radial to the tendon of the flexor carpi radialis (FCR). Incise forearm fascia between FCR and vascular bundle of nerves. Go to the pronator quadratus. Make an L-shaped incision on the fascia of the pronator quadratus and push the muscle from the distal radius towards the ulna. Then set the fracture.

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
POLYAXIAL IMPLANTATION

Polyaxial implantation of the Polyaxial Winkelstabil distal radius plate
Example with right side plate and polyaxial screw-in drill sleeve

**POLYAXIAL WINKELSTABIL SINGLE USE INSTRUMENTE UND IMPLANTATE
POLYAXIAL ANGLE STABLE SINGLE USE INSTRUMENTS AND IMPLANTS**

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POLYAXIALE IMPLANTATION

6. Inserting the angle-stable screws

Always start with the distal row. The selection of the first hole depends on the progression direction of the fracture. Screw the polyaxial drill sleeve into the desired plate hole. The funnel of the drill sleeve allows for a free angle of up to 15°. Check the direction of the drill with the image converter to avoid damage to the joint surface from the drill.

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
POLYAXIAL IMPLANTATION

12. Rear distal row of screws

If necessary, additional screws can be inserted into the rear distal screw hole plane. Inserting the screws into the rear distal screw hole plane does not necessarily require an image converter check during drilling.

13. Filling the proximal shaft holes

When the distal fragments are sufficiently secured, the next step is to finish fastening the plate in the proximal shaft section. To do so, screw the drill sleeve into the desired proximal plate hole and carefully drill the hole through both corticals.

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

**MONOAXIAL IMPLANTATION WITH SCREW-IN
DRILL SLEEVE**

Monoaxial implantation of the Polyaxial Winkelstabil
distal radius plate
Example with right side plate and monoaxial drill sleeve

**POLYAXIAL WINKELSTABIL SINGLE USE INSTRUMENTE UND IMPLANTATE
POLYAXIAL ANGLE STABLE SINGLE USE INSTRUMENTS AND IMPLANTS**

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1. Temporary fixation of the plate

After the reduction of the fracture, place the plate on the volar radius surface and temporarily affix plate by threading the 1.6 mm Kirschner wires through the Kirschner wire holes in the bone. Check the correct positioning of the plate by means of an image converter.

2. Filling the oblong hole

Position the drill sleeve in the center of the oblong hole and drill carefully through both corticals using the 1.9 mm drill.

3. Measuring the screw length

Use the depth gauge to determine the required screw length.

4. Screwing in of the TX screw

Now insert the TX screw into the plate hole. Do not completely tighten the screw yet.

5. Removing the K-wire

Now remove the K-wire. The oblong hole makes it possible to adjust the plate slightly to achieve the correct position. Check the correct plate position with the image converter. Once the plate is in the correct position, it is now advisable to tighten the TX screw for fixing the plate position in place.

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MONOAXIAL IMPLANTATION WITH SCREW-IN DRILL SLEEVE

7. Determining the screw length

Remove the drill sleeve and measure the required screw length with the depth measuring gauge.

8. Screwing in the PTXL screw

Now screw the PTXL screw into the plate hole. The screw heads of the PTXL screws are conical in shape and must therefore only be fastened hand-tight. Check the correct position of the screw with the image converter.

9. Filling the additional distal holes

Screw the drill sleeve into the next plate hole to be filled. Use the image converter to check the direction of the drill to avoid damage to the joint surface.

10. Determining the screw length

Remove the drill sleeve and measure the required screw length using the depth gauge.

11. Screwing in additional PTXL screws

Now screw the PTXL screw into the plate hole. The screw heads of the PTXL screws are conical in shape and must therefore only be fastened hand-tight. Check the correct position of the screw with the image converter.

12. Rear distal row of screws

If necessary, additional screws can be inserted into the rear distal screw hole plane.

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MONOAXIAL IMPLANTATION WITH SCREW-IN DRILL SLEEVE

13. Filling the proximal shaft holes

When the distal fragments are sufficiently secured, the next step is to finish fastening the plate in the proximal shaft section. To do so, screw the drill sleeve into the desired proximal hole and carefully drill the hole through both corticals.

14. Determining the screw length

Determine the screw length same as in the previously described steps.

15. Screwing the PTXL screw in the shaft

Except for the oblong hole, all other proximal screw holes are filled with angle-stable PTXL screws. It is not necessary to fill all proximal holes with screws.

16. Completion

Once you have inserted all the necessary screws for the fixation of the fracture, it is best to check the result by taking a final image with the image converter.

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RETRAX, YOUR CLEVER SUPPORT FOR ALL RADIUS FRACTURES

Self-retaining and self-stabilizing single-use retractor system for the treatment of radius fractures

Retrax offers many advantages:

- » Easy-to-use retractor system offering a very good view of the surgical site.
- » Two forearm braces to ensure the optimal forearm fit.
- » All retractors and braces are made from recyclable synthetic material.
- » The synthetic material is extremely stable, but also highly flexible.
- » The retractor remains in place during the entire surgical procedure and offers extremely stable access.
- » Retrax is 100% x-ray permeable.
- » Retrax eliminates the need to remove metal hooks for image converter check during the operation.
- » High stability allows you to keep your forearm far from the x-ray beam cone. This means that you reduce your own direct x-ray exposure.
- » The system enables you to treat radius fractures even without an assistant, in an emergency.
- » If used correctly, Retrax helps to reduce surgery time.


Contraindications:

- » Patients who are very thin or strongly overweight, on whom the arm brace does not fit correctly.
- » Skin conditions that make it vulnerable (e.g. parchment skin from cortisone use).
- » Vascular diseases causing reduced blood supply (circulation around the surgical site is disturbed).
- » Excessive duration of surgery (150 min.).
- » Any medical condition or surgical situation posing a contraindication to the use of the retractor system.

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
HOW TO USE RETRAX

1. Select the appropriate arm brace from the kit and place the forearm in the brace. In doing so, please pay attention to the position of the „distal“ label.
2. The proximal fixation element should preferably be inserted into the most proximal fixation slot of the arm brace and then locked using the click mechanism.
3. Cut into the skin. Open and spread the subcutaneous tissue.
4. Retract the subcutaneous tissue with one or two soft-tissue hooks.
5. Prepare the tissue to the bone and create soft tissue canals for the tips of the Hohmann retractor at the bone margin with dissecting scissors.
6. Insert the Hohmann retractors at the chosen sites and lock them into position, too, using the click mechanism.

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




HOW TO USE RETRAX

7. Insert a soft tissue hook through a slot on the proximal fixation element and retract it proximally.
8. Insert the fixing clamp into the chosen interdigital space.
9. Insert another soft tissue hook into the slot on the fixation clamp and position it in the wound opening. Retract it distally.
10. Reposition and place a plate. The plate holder is used to fix it. Carry out an x-ray check with Retrax in a horizontal position and with the temporary plate holder and provisionally fixed plate.
11. Lock the correctly aligned plate into place with the first screw.
12. Remove the plate holder and place screws in the remaining holes. Final x-ray check. Remove Retrax and dispose of it in an environmentally appropriate manner.

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POLYAXIAL WINKELSTABIL FIBULA PLATTEN LINKS


Polyaxial angle stable DISTAL FIBULA PLATE

Perfect design and 100% made in Germany

Design of the fibula plate

The anatomic design of the Polyaxial Winkelstabil fibula plate is unique. The tip on the distal end of the plate can be anchored very well in osteoporotic bone. If the tip is not needed, it can be pinched off with a side cutter. Distally up to seven screws can be screwed in in a angle-stable fashion in a balanced direction. Fibula that have been fractured into small pieces are reconstructed in the spoon-like design of the plate. The high degree of stability enables rapid weight-bearing mobilization in comparison to conventional systems.

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Polyaxial Winkelstabil DISTAL FIBULA PLATE
Perfect design and 100% made in Germany

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Indication

The 2.7 mm angle-stable distal fibula system is indicated for the treatment of isolated fibula or fibula fractures in the context of complex upper ankle joint injuries. The various screw connection options in the distal portion offer advantages especially in case of severe osteoporosis.

Monoaxial stabilization

The anatomical design of the Polyaxial angle stable fibula plate is unique. The angles of the screws pre-set in the plate are perfectly adapted to the anatomical conditions of the fibula. Therefore, we recommend the implantation of angle-stable PTXL screws only in conjunction with the monoaxial screw-in drill sleeve.

Distal mandrel of the fibula plate

The mandrel at the distal end of the plate can be well anchored in osteoporotic bone. If the tip is not required, it can be pinched off with a wire cutter.

Self-locking screws

The Polyaxial angle stable T10 screwdriver is manufactured with very high precision. This high precision allows for the screwdriver tip being „wedged“ in the screw head so that it is firmly connected with the screwdriver. When positioning the screwdriver in the head of the screw, make sure to press the screwdriver down into the screw head using a certain amount of pressure. If using insufficient pressure, the tip of the screwdriver may not penetrate deep enough into the screw head. If this is the case, the screw head may become damaged during tightening. Polyaxial angle stable more than just plates and screws

Sterile implants

The Polyaxial angle stable implants are all packaged individually and sterile and are color-coded. Please note the correct side when selecting the plate. Check packaging and labeling on the plate to ensure that the plate is for the correct side („R“ for right and „L“ for left).

Soft tissue procedure

A longitudinal incision across the fibula, starting at the outer ankle top, about 10-12 cm in length towards the proximal end, provides a good overview of the site for an anatomical reduction and placement of the plate. In the distal half of the incision, it is possible to cut directly on the fibula. A layered preparation is advisable for the proximal sections to protect the variably progressing portion of the peroneal nerve.

Positioning the plate

After reduction of the fracture, position the plate distally with the mandrel first. If the mandrel proves to be disruptive, it can be pinched off with pliers. Fit a short TX screw into the most distal hole but do not finish tightening it yet. Fine tune the placement of the plate around this fixation so that the most proximal hole is securely positioned over the fibula. If needed, a Kirschner wire can be placed in the most proximal hole. With this, the plate is already relatively reliably fixed in place.

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Polyaxial Winkelstabil single use Instrumente und
Implantate

FIBULA PLATE IMPLANTATION

Monoaxial implantation of the Polyaxial Winkelstabil
fibula plate
Example with left side plate and screw-in drill sleeve

Polyaxial angle stable single use instruments and
implants

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IMPLANTATION FIBULA PLATE

1. Temporary fixation of the plate, distal
After the reduction of the fracture, first press the mandrel into the distal end of the fibula.
2. Temporary fixation of the plate, proximal
Then position the plate at the proximal end and fix this position with a K-wire.
3. Distal screw stabilization
As a rule, it is advisable to first secure the distal fragments. To do so, insert the drill sleeve into the desired plate hole and drill the hole for the PTXL screw above it.
4. Measuring the screw length
Use the measuring gauge to measure the screw length.
5. Screwing in the PTXL screw
Now screw the PTXL screw into the plate hole. The screw heads of the PTXL screws are conical in shape and must therefore only be fastened hand-tight. Check the correct position of the screw with the image converter.
6. Screwing in additional PTXL screws
Repeat the procedure for additional PTXL screws you wish to place in the distal area until the distal fragments are sufficiently fixed.

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IMPLANTATION FIBULA PLATE

7. PTXL screws in the proximal area

After stabilization of the distal fragments, the next step consists of fixing the plate on the proximal end. To do so, screw the drill sleeve into the desired plate hole. Drilling and measuring is carried out the same as described previously.

8. Verification of the syndesmosis ligament

Check the integrity of the syndesmosis either clinically or by means of the image converter. In case of instability, it is possible to fit a screw. The long TX screws are used to fit the adjustment screw.

If the bone quality is very poor, a variant that can be used in exceptional cases consists of a multiple screw connection of the fibula against the tibia by also using the ventral PTXL screws of the plate, proximally of the oblong hole. Please note that this advanced procedure is not considered standard treatment and that it is recommended to separately justify and explain the reason for using this type of treatment in the surgical report.

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Polyaxial angle stable SINGLE USE RADIUS SYSTEM*

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Polyaxial angle stable SINGLE USE FIBULA SYSTEM*
Basic recommended items

Art.Nr. / Ref.No.	Beschreibung / Description	Menge / Quantity
93-7100-020H	single use R+F Instrumenten Kit f. 2.7 mm T10 Schrauben single use R+F instruments kit f. 2.7 mm, T10 screws	1
93-6600-016	DF-Platte 7/6 Loch rechts DF-plate 7/6 hole right	1
93-6600-017	DF-Platte 7/6 Loch links DF-plates 7/6 hole left	1

93-6127-010TX	Kort. Tit. TX-Schr. Ø 2.7 mm, L. 10 mm, T10 Cort. tit. TX-screw, Ø 2.7 mm, L. 10 mm, T10	1
93-6127-012TX	Kort. Tit. TX-Schr. Ø 2.7 mm, L. 12 mm, T10 Cort. tit. TX-screw, Ø 2.7 mm, L. 12 mm, T10	1
93-6127-014TX	Kort. Tit. TX-Schr. Ø 2.7 mm, L. 14 mm, T10 Cort. tit. TX-screw, Ø 2.7 mm, L. 14 mm, T10	1
93-6127-016TX	Kort. Tit. TX-Schr. Ø 2.7 mm, L. 16 mm, T10 Cort. tit. TX-screw, Ø 2.7 mm, L. 16 mm, T10	1
93-6127-038TX	Kort. Tit. TX-Schr. Ø 2.7 mm, L. 38 mm, T10 Cort. tit. TX-screw, Ø 2.7 mm, L. 38 mm, T10	1
93-6127-040TX	Kort. Tit. TX-Schr. Ø 2.7 mm, L. 40 mm, T10 Cort. tit. TX-screw, Ø 2.7 mm, L. 40 mm, T10	1
93-6127-044TX	Kort. Tit. TX-Schr. Ø 2.7 mm, L. 44 mm, T10 Cort. tit. TX-screw, Ø 2.7 mm, L. 44 mm, T10	1

93-6127-010PTXL	Kort. Tit. PTXL-Schraube, Ø 2.7 mm, L. 10 mm, T10 Cort. tit. PTXL-screw, Ø 2.7 mm, L. 10 mm, T10	2
93-6127-012PTXL	Kort. Tit. PTXL-Schraube, Ø 2.7 mm, L. 12 mm, T10 Cort. tit. PTXL-screw, Ø 2.7 mm, L. 12 mm, T10	4
93-6127-014PTXL	Kort. Tit. PTXL-Schraube, Ø 2.7 mm, L. 14 mm, T10 Cort. tit. PTXL-screw, Ø 2.7 mm, L. 14 mm, T10	4
93-6127-016PTXL	Kort. Tit. PTXL-Schraube, Ø 2.7 mm, L. 16 mm, T10 Cort. tit. PTXL-screw, Ø 2.7 mm, L. 16 mm, T10	3
93-6127-018PTXL	Kort. Tit. PTXL-Schraube, Ø 2.7 mm, L. 18 mm, T10 Cort. tit. PTXL-screw, Ø 2.7 mm, L. 18 mm, T10	3
93-6127-020PTXL	Kort. Tit. PTXL-Schraube, Ø 2.7 mm, L. 20 mm, T10 Cort. tit. PTXL-screw, Ø 2.7 mm, L. 20 mm, T10	2
93-6127-022PTXL	Kort. Tit. PTXL-Schraube, Ø 2.7 mm, L. 22 mm, T10 Cort. tit. PTXL-screw, Ø 2.7 mm, L. 22 mm, T10	2
93-6127-024PTXL	Kort. Tit. PTXL-Schraube, Ø 2.7 mm, L. 24 mm, T10 Cort. tit. PTXL-screw, Ø 2.7 mm, L. 24 mm, T10	2

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