



STERIL CONTAINER SYSTEM + SIEBKÖRBE

STERILE CONTAINER SYSTEMS + BASKETS



Einleitung

Die Sterilcontainer Systeme ermöglichen eine systematische Organisation Ihres gesamten Sterilgutes- Die komfortable Handhabung und die Vielfalt unserer Container bietet Ihnen Sicherheit und Kontrolle in allen Anwendungsbereichen der Sterilgutversorgung- Lagerung, Transport und Entsorgung können mit unseren Containern perfekt organisiert werden- Unsere Sterilcontainer entsprechen den Anforderungen der Normen und Standards zur Sterilgutverpackung- Die Container sind hergestellt aus eloxiertem Aluminium und zeichnen sich somit durch geringes Gewicht und gutes Trocknungsverhalten aus- Farbige Deckel und Kennzeichnungsschilder ermöglichen Ihnen eine genaue Kontrolle über Ihr Sterilgut- Durch unser umfangreiches Programm bieten wir Ihnen eine genaue Anpassung an Ihre Bedürfnisse-

Introduction

The sterilisation container systems allow you a systematic organisation of your entire instrument sterilisation- The comfortable handling and diversity of our containers gives you safety and control for your instrument sterilisation- Storage, transportation and disposal are perfectly organised with our containers- Our sterilisation container systems are conform to the international standard requirements- The sterilisation containers are made of anodised aluminium, have minimum weight and excellent drying properties- Colour-coded lids and identification labels enable you precise controlling of your instrument sterilisation- Any changes or adjustments in your requirements are no problem due to our extensive range of products-

Suitability of our sterile containers for steam sterilization procedures¹⁾ as per DIN 58953, part 9 (issue May 1987)

Sterilization container	Packing		Steam sterilization procedure acc.to DIN 58946, part I				
	Description	Symbolic representation and packing system	Gravitation procedure	Pre - vacuum procedure	Fractional vacuum procedure	Fractional flow procedure	Steam injection procedure
with filters (StB - F)	Textiles in sterilization container with filter in lid and bottom acc.to DIN 58952 part I		+	+	+	+	+
			0	0	+	+	+
	Textiles in sterilization container with filter in lid, bottom not perforated acc.to DIN 58952 part I		—	—	+	+	+
			—	—	+	+	0
	Textiles in sterilization container with filter in bottom lid, not perforated acc.to DIN 58952 part I		0	+	+	+	+
			—	0	+	+	+
	Instruments in wire basket acc.to DIN 58952 part III in sterilization container with filter in lid and bottom acc. to DIN 58952 part I		+ ³⁾	+	+	+	+
	Instruments in wire basket acc.to DIN 58952 part III, in sterilization container with filter in lid and bottom not perforated acc. to DIN 58952 part I		—	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket acc.to DIN 58952 part III, in sterilization container with filter in bottom acc. to DIN 58952 part I		0	+	+	+	+
with valves (StB - V)	Textiles in sterilization container with valves in lid bottom acc. to DIN 58952 part I		—	0	+	+	0
			—	—	+	+	0
	Textiles in sterilization container with valves in lid no valves on the bottom		—	—	+	+	—
			—	—	+	+	—
	Instruments in wire basket in sterilization container with valve in lid and bottom		—	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾
	Instruments in wire basket in sterilization container with valve in lid and not perforated		—	+ ²⁾	+ ²⁾	+ ²⁾	+ ²⁾

Explanation of symbols : + applicable — not applicable 0 limited application

1) In order to avoid mistakes with regards to service, operation and/or when loading the sterilizer, please refer to table 1 for correct or wrong combination of sterilization containers for standardized steam sterilization procedures.
2) Only sterilization containers guaranteeing evaporation or discharge of condensate are suitable above.
3) Gravitation procedures having no drying device do not guarantee a sufficient after - drying effect.

Inhaltsverzeichnis / Index

Seite / Page

3	Einleitung / Preface
4	Tabelle / Table
6 - 7	Information / Information
9 - 13	1/1 Container und Siebkörbe / 1/1 Container and Wire Baskets
15 - 19	3/4 Container und Siebkörbe / 3/4 Container and Wire Baskets
21 - 25	1/2 Container und Siebkörbe / 1/2 Container and Wire Baskets
27 - 29	Großraumcontainer und Siebkörbe / Wide Body Container and Wire Baskets
31 - 33	Flat Container und Siebkörbe / Flat Container and Wire Baskets
34 - 35	Kleinsetcontainer / Mini Container
36	Dental Container
37 - 38	Implantat und Optikcontainer / Implant and Scopy Container
39 - 40	Entsorgungscontainer / Removal Container
41	Bürstencontainer / Brush Container
43 - 46	Zubehör / Accessories
47	Dauerfilter / Permanent Filter
48-50	Edelstahlsiebkörbe / Trays made of stainless steel
51	Silikonlagerungen / Silicon fixation
52 - 53	Side Perforated Baskets Wire Base / Lid for Side Perforated Baskets Wire Base
54 - 55	Side Perforated Basket Flat Base / Lid for Side Perforated Basket Flat Base
56 - 57	Standard Wire Mesh Baskets / Lid for Standard Wire Mesh Baskets
58 - 59	Crimped Wire Mesh Baskets / Lids for Crimped Wire Mesh Baskets
60 - 61	Classic Crimped Wire Mesh Baskets / Fine Mesh Baskets With Lid
62	Micro Fine Mesh Baskets With Lid
63	Fine Mesh Baskets with Unique Twist-Off Lid
64	Micro Fine Mesh Baskets with Unique Twist-Off Lid
65	Full Perforated Trays with Lid
66	Wire Net Baskets
67	Sterile Goods Baskets
68	BASKET SPRI
69	Wire Racks
70 - 72	Endoscope Baskets
73	Micro Ophthalmic Baskets
74	Micro Instruments Cassettes
75 - 76	Instruments Cassettes
77	Micro Cassettes Fixing Clamps
78 - 79	Implantology Cassettes
80	Bur Clips / Bur Holder Stand
81 - 83	Hospital Ware
84	Hospital Ware / Basket Accessories / Silicon Matt & Holders
85	Notizen



Wählen Sie zwischen unserem Standard-, und unserem Sicherheitsmodell-
Der Sicherheitsdeckel bietet zusätzlichen Schutz vor Kontamination bei Transport und Lagerung- Filter-, und Sicherheitsdeckel lassen sich nur bei geöffnetem Behälter trennen-

Choose between our standard and our model-
The safety lid protects against contamination during transportation and storage- The filter- and safety-lid can only be separated by complete opening of the container-





Die Verschlusslaschen sind in jeder Position arretierbar und bis 180° in den Deckel einschwenkbar-

The locking flaps are movable 180° on a pivot and can be secured in any position-



Die Filterhalterung ist mit einfachem Druck auf die Filtermitte einzusetzen, mit Druck auf die Lasche lösen Sie die Filterhalterung wieder-

The filter retainer is locked by pressing down in center of the filter- Opening the filter retainer by pressing the plate

1/1 Container 1/1 Size Containers



Boden ungelocht
Bottom non perforated



	99-110-100
	99-120-100
	99-130-100
	99-140-100
	99-150-100
	99-160-100

Gewicht · weight
2634 g



außen · outside	580 x 280 x 100 mm
innen · inside	550 x 265 x 84 mm

Boden gelocht
Bottom perforated



	99-111-100
	99-121-100
	99-131-100
	99-141-100
	99-151-100
	99-161-100

Gewicht · weight
3156 g

	99-110-135
	99-120-135
	99-130-135
	99-140-135
	99-150-135
	99-160-135

Gewicht · weight
3129,5 g



außen · outside	580 x 280 x 135 mm
innen · inside	550 x 265 x 115 mm

	99-111-135
	99-121-135
	99-131-135
	99-141-135
	99-151-135
	99-161-135

Gewicht · weight
3244 g

	99-110-150
	99-120-150
	99-130-150
	99-140-150
	99-150-150
	99-160-150

Gewicht · weight
3219 g



außen · outside	580 x 280 x 150 mm
innen · inside	550 x 265 x 130 mm

	99-111-150
	99-121-150
	99-131-150
	99-141-150
	99-151-150
	99-161-150

Gewicht · weight
3382 g

	99-110-200
	99-120-200
	99-130-200
	99-140-200
	99-150-200
	99-160-200

Gewicht · weight
3429,5 g



außen · outside	580 x 280 x 200 mm
innen · inside	550 x 265 x 180 mm

	99-111-200
	99-121-200
	99-131-200
	99-141-200
	99-151-200
	99-161-200

Gewicht · weight
3738,5 g

	99-110-260
	99-120-260
	99-130-260
	99-140-260
	99-150-260
	99-160-260

Gewicht · weight
3919 g



außen · outside	580 x 280 x 260 mm
innen · inside	550 x 265 x 240 mm

	99-111-260
	99-121-260
	99-131-260
	99-141-260
	99-151-260
	99-161-260

Gewicht · weight
4272,5 g

1/1 Size Container

Safe Model · Sicherheitsmodell

Boden ungelocht
Bottom non perforated




	99-112-100
	99-122-100
	99-132-100
	99-142-100
	99-152-100
	99-162-100

Gewicht · weight
3486 g



außen · outside	580 x 280 x 100 mm
innen · inside	550 x 265 x 84 mm

Boden gelocht
Bottom perforated



	99-113-100
	99-123-100
	99-133-100
	99-143-100
	99-153-100
	99-163-100

Gewicht · weight
4007,5 g

	99-112-135
	99-122-135
	99-132-135
	99-142-135
	99-152-135
	99-162-135

Gewicht · weight
3981,5 g



außen · outside	580 x 280 x 135 mm
innen · inside	550 x 265 x 115 mm

	99-113-135
	99-123-135
	99-133-135
	99-143-135
	99-153-135
	99-163-135

Gewicht · weight
4095,5 g

	99-112-150
	99-122-150
	99-132-150
	99-142-150
	99-152-150
	99-162-150

Gewicht · weight
4070,5 g



außen · outside	580 x 280 x 150 mm
innen · inside	550 x 265 x 130 mm

	99-113-150
	99-123-150
	99-133-150
	99-143-150
	99-153-150
	99-163-150

Gewicht · weight
4233 g

	99-112-200
	99-122-200
	99-132-200
	99-142-200
	99-152-200
	99-162-200

Gewicht · weight
4280 g



außen · outside	580 x 280 x 200 mm
innen · inside	550 x 265 x 180 mm

	99-113-200
	99-123-200
	99-133-200
	99-143-200
	99-153-200
	99-163-200

Gewicht · weight
4589 g

	99-112-260
	99-122-260
	99-132-260
	99-142-260
	99-152-260
	99-162-260

Gewicht · weight
4770,5 g



außen · outside	580 x 280 x 260 mm
innen · inside	550 x 265 x 240 mm


	99-113-260
	99-123-260
	99-133-260
	99-143-260
	99-153-260
	99-163-260

Gewicht · weight
5140 g

Boden ungelocht
Bottom non perforated




Boden gelocht
Bottom perforated



	außen · outside	innen · inside	
99-105-100	580 x 280 x 100 mm	550 x 262 x 84 mm	1923 g
99-105-135	580 x 280 x 135 mm	550 x 262 x 117 mm	2019 g
99-105-150	580 x 280 x 150 mm	550 x 262 x 132 mm	2108 g
99-105-200	580 x 280 x 200 mm	550 x 262 x 182 mm	2317 g
99-105-260	580 x 280 x 260 mm	550 x 262 x 244 mm	2808 g

	außen · outside	innen · inside	
99-106-100	580 x 280 x 100 mm	550 x 262 x 84 mm	2044 g
99-106-135	580 x 280 x 135 mm	550 x 262 x 117 mm	2132 g
99-106-150	580 x 280 x 150 mm	550 x 262 x 132 mm	2270 g
99-106-200	580 x 280 x 200 mm	550 x 262 x 182 mm	2627 g
99-106-260	580 x 280 x 260 mm	550 x 262 x 244 mm	2851 g

Deckel gelocht
Lid perforated



580 x 280 mm

Sicherheitsdeckel
Safety Lid




	99-100-010
	99-100-020
	99-100-030
	99-100-040
	99-100-050
	99-100-060


Gewicht · weight
1111 g

	99-103-010
	99-103-020
	99-103-030
	99-103-040
	99-103-050
	99-103-060


Gewicht · weight
851,5 g



Siebkorb · Wire Basket		
99-190-030	485 x 255 x 30 mm	1108 g
99-190-050	485 x 255 x 50 mm	1112 g
99-190-070	485 x 255 x 70 mm	1176 g
99-190-100	485 x 255 x 100 mm	1289 g



Siebkorb · Wirebasket		
99-191-030	540 x 255 x 30 mm	1149 g
99-191-050	540 x 255 x 50 mm	1189 g
99-191-070	540 x 255 x 70 mm	1259 g
99-191-100	540 x 255 x 100 mm	1354 g



Zwischenboden · Wire Base		
99-193-0000	545 x 255 x 30 mm	729 g



Spannrahmen · Drape Retainer		
99-195-000	545 x 255 mm	




Einschlagtuch · Drape	
99-095-070	130 x 60 cm
99-095-080	140 x 100 cm
99-095-090	185 x 45 cm
99-095-100	185 x 60 cm

3/4 Container 3/4 Size Containers



Boden ungelocht
Bottom non perforated



	99-210-100
	99-220-100
	99-230-100
	99-240-100
	99-250-100
	99-260-100

Gewicht · weight
2268,5 g



außen · outside	465 x 280 x 100 mm
innen · inside	425 x 260 x 84 mm

Boden gelocht
Bottom perforated



	99-211-100
	99-221-100
	99-231-100
	99-241-100
	99-251-100
	99-261-100

Gewicht · weight
2503 g

	99-210-135
	99-220-135
	99-230-135
	99-240-135
	99-250-135
	99-260-135

Gewicht · weight
2549,5 g



außen · outside	465 x 280 x 135 mm
innen · inside	425 x 260 x 115 mm

	99-211-135
	99-221-135
	99-231-135
	99-241-135
	99-251-135
	99-261-135

Gewicht · weight
2715 g

	99-210-150
	99-220-150
	99-230-150
	99-240-150
	99-250-150
	99-260-150

Gewicht · weight
2620,5 g



außen · outside	465 x 280 x 150 mm
innen · inside	425 x 260 x 130 mm

	99-211-150
	99-221-150
	99-231-150
	99-241-150
	99-251-150
	99-261-150

Gewicht · weight
2718 g

3/4 Size Container

Sicherheitsmodell · Safe Model

Boden ungelocht
Bottom non perforated



	99-212-100
	99-222-100
	99-232-100
	99-242-100
	99-252-100
	99-262-100

Gewicht · weight
2980 g



außen · outside	465 x 280 x 100 mm
innen · inside	425 x 260 x 84 mm

Boden gelocht
Bottom perforated



	99-213-100
	99-223-100
	99-233-100
	99-243-100
	99-253-100
	99-263-100

Gewicht · weight
3213,5 g

	99-212-135
	99-222-135
	99-232-135
	99-242-135
	99-252-135
	99-262-135

Gewicht · weight
3260 g



außen · outside	465 x 280 x 135 mm
innen · inside	425 x 260 x 115 mm

	99-213-135
	99-223-135
	99-233-135
	99-243-135
	99-253-135
	99-263-135

Gewicht · weight
2425 g

	99-212-150
	99-222-150
	99-232-150
	99-242-150
	99-252-150
	99-262-150

Gewicht · weight
2975 g



außen · outside	465 x 280 x 150 mm
innen · insid	425 x 260 x 130 mm


	99-213-150
	99-223-150
	99-233-150
	99-243-150
	99-253-150
	99-263-150

Gewicht · weight
3561 g

Boden ungelocht
Bottom non perforated




Boden gelocht
Bottom perforated



	außen · outside	innen · inside	
99-205-100	465 x 280 x 100 mm	425 x 262 x 85 mm	1396 g
99-205-135	465 x 280 x 135 mm	425 x 262 x 121 mm	1676 g
99-205-150	465 x 280 x 150 mm	425 x 262 x 135 mm	1792 g

	außen · outside	innen · inside	
99-206-100	465 x 280 x 100 mm	425 x 262 x 85 mm	1629 g
99-206-135	465 x 280 x 135 mm	425 x 262 x 121 mm	1842 g
99-206-150	465 x 280 x 150 mm	425 x 262 x 135 mm	1978 g

Deckel gelocht
Lid perforated



465 x 280 mm

Sicherheitsdeckel
Safety Lid




	99-200-010
	99-200-020
	99-200-030
	99-200-040
	99-200-050
	99-200-060


Gewicht · weight
873 g

	99-203-010
	99-203-020
	99-203-030
	99-203-040
	99-203-050
	99-203-060

Gewicht · weight
710 g



Drahtsiebkorb · Wire Basket		
99-290-030	405 x 255 x 30 mm	972 g
99-290-050	405 x 255 x 50 mm	988 g
99-290-070	405 x 255 x 70 mm	1061 g
99-290-100	405 x 255 x 100 mm	1143 g



Zwischenboden · Wire Base		
99-293-030	405 x 255 x 30 mm	598 g



Spannrahmen · Drape Retainer	
99-295-000	405 x 255 mm




Einschlagtuch · Drape	
99-095-080	140 x 100 cm

1/2 Container 1/2 Size Containers



Boden ungelocht
Bottom non perforated



	99-310-100
	99-320-100
	99-330-100
	99-340-100
	99-350-100
	99-360-100

Gewicht · weight
1744,5 g



außen · outside	285 x 280 x 100 mm
innen · inside	262 x 260 x 84 mm

Boden gelocht
Bottom perforated



	99-311-100
	99-321-100
	99-331-100
	99-341-100
	99-351-100
	99-361-100

Gewicht · weight
1890,5 g

	99-310-135
	99-320-135
	99-330-135
	99-340-135
	99-350-135
	99-360-135

Gewicht · weight
1947 g



außen · outside	285 x 280 x 135 mm
innen · inside	262 x 260 x 116 mm

	99-311-135
	99-321-135
	99-331-135
	99-341-135
	99-351-135
	99-361-135

Gewicht · weight
2125 g

	99-310-150
	99-320-150
	99-330-150
	99-340-150
	99-350-150
	99-360-150

Gewicht · weight
1981 g



außen · outside	285 x 280 x 150 mm
innen · inside	262 x 260 x 132 mm

	99-311-150
	99-321-150
	99-331-150
	99-341-150
	99-351-150
	99-361-150

Gewicht · weight
2462,5 g

	99-310-200
	99-320-200
	99-330-200
	99-340-200
	99-350-200
	99-360-200

Gewicht · weight
2101 g



außen · outside	285 x 280 x 200 mm
innen · inside	262 x 260 x 182 mm

	99-311-200
	99-321-200
	99-331-200
	99-341-200
	99-351-200
	99-361-200

Gewicht · weight
2513,5 g

	99-310-260
	99-320-260
	99-330-260
	99-340-260
	99-350-260
	99-360-260

Gewicht · weight
2483,5 g



außen · outside	285 x 280 x 260 mm
innen · inside	262 x 260 x 238 mm

	99-311-260
	99-321-260
	99-331-260
	99-341-260
	99-351-260
	99-361-260

Gewicht · weight
2657,5 g

1/2 Size Container

Sicherheitsmodell · Safe Model

Boden ungelocht
Bottom non perforated



	99-312-100
	99-322-100
	99-332-100
	99-342-100
	99-352-100
	99-362-100

Gewicht · weight
2162,5 g



außen · outside	285 x 280 x 100 mm
innen · inside	262 x 260 x 84 mm

Boden gelocht
Bottom perforated



	99-313-100
	99-323-100
	99-333-100
	99-343-100
	99-353-100
	99-363-100

Gewicht · weight
2309,5 g

	99-312-135
	99-322-135
	99-332-135
	99-342-135
	99-352-135
	99-362-135

Gewicht · weight
2351 g



außen · outside	285 x 280 x 135 mm
innen · inside	262 x 260 x 116 mm

	99-313-135
	99-323-135
	99-333-135
	99-343-135
	99-353-135
	99-363-135

Gewicht · weight
2466 g

	99-312-150
	99-322-150
	99-332-150
	99-342-150
	99-352-150
	99-362-150

Gewicht · weight
2394 g



außen · outside	285 x 280 x 150 mm
innen · inside	262 x 260 x 132 mm

	99-313-150
	99-323-150
	99-333-150
	99-343-150
	99-353-150
	99-363-150

Gewicht · weight
2583 g

	99-312-200
	99-322-200
	99-332-200
	99-342-200
	99-352-200
	99-362-200

Gewicht · weight
2560 g



außen · outside	285 x 280 x 200 mm
innen · inside	262 x 260 x 182 mm

	99-313-200
	99-323-200
	99-333-200
	99-343-200
	99-353-200
	99-363-200

Gewicht · weight
2785,5 g

	99-312-260
	99-322-260
	99-332-260
	99-342-260
	99-352-260
	99-362-260

Gewicht · weight
2930 g



außen · outside	285 x 280 x 260 mm
innen · inside	262 x 260 x 238 mm


	99-313-260
	99-323-260
	99-333-260
	99-343-260
	99-353-260
	99-363-260

Gewicht · weight
3069,5 g

Boden ungelocht
Bottom non perforated




Boden gelocht
Bottom perforated



	außen · outside	innen · inside	
99-305-100	285 x 280 x 100 mm	262 x 262 x 82 mm	1096 g
99-305-135	285 x 280 x 135 mm	262 x 262 x 120 mm	1284 g
99-305-150	285 x 280 x 150 mm	262 x 262 x 135 mm	1327 g
99-305-200	285 x 280 x 200 mm	262 x 262 x 185 mm	1447 g
99-305-260	285 x 280 x 260 mm	262 x 262 x 240 mm	1817 g

	außen · outside	innen · inside	
99-306-100	285 x 280 x 100 mm	262 x 262 x 82 mm	1242 g
99-306-135	285 x 280 x 135 mm	262 x 262 x 120 mm	1398 g
99-306-150	285 x 280 x 150 mm	262 x 262 x 135 mm	1471 g
99-306-200	285 x 280 x 200 mm	262 x 262 x 185 mm	1672 g
99-306-260	285 x 280 x 260 mm	262 x 262 x 240 mm	2000 g

Deckel gelocht
Lid perforated



285 x 280 mm

Sicherheitsdeckel
Safety Lid



	99-300-010
	99-300-020
	99-300-030
	99-300-040
	99-300-050
	99-300-060


Gewicht · weight
654 g

	99-303-010
	99-303-020
	99-303-030
	99-303-040
	99-303-050
	99-303-060

Gewicht · weight
458 g



Drahtsiebkorb · Wire Basket		
99-390-030	255 x 245 x 30 mm	637,5 g
99-390-050	255 x 245 x 50 mm	648,5 g
99-390-070	255 x 245 x 70 mm	691,5 g
99-390-100	255 x 245 x 100 mm	758,5 g



Zwischenboden · Wire Base		
99-393-000	258 x 258 x 30 mm	351 g



Spannrahmen · Drape Retainer	
99-395-000	255 x 255 mm




Einschlagtuch · Drape	
99-095-040	90 x 60 cm
99-095-050	90 x 80 cm
99-095-060	130 x 45 cm
99-095-070	130 x 60 cm

Großraum Container Wide Body Containers



Boden ungelocht
Bottom non perforated



	99-110-62
	99-120-62
	99-130-62
	99-140-62
	99-150-62
	99-160-62

Gewicht · weight
3857 g



outside	600 x 400 x 120 mm
inside	550 x 350 x 100 mm

Boden gelocht
Bottom perforated



	99-111-62
	99-121-62
	99-131-62
	99-141-62
	99-151-62
	99-161-62

Gewicht · weight
4101 g




outside	600 x 400 x 180 mm
inside	550 x 350 x 160 mm

	99-110-68
	99-120-68
	99-130-68
	99-140-68
	99-150-68
	99-160-68


Gewicht · weight
4043,5 g

	99-111-68
	99-121-68
	99-131-68
	99-141-68
	99-151-68
	99-161-68

Gewicht · weight
4581,5 g



Drahtsiebkorb · Wire Basket	
99-190-55	470 x 337 x 51 mm



Drahtsiebkorb · Wire Basket	
99-191-60	533 x 337 x 101 mm

Boden ungelocht
Bottom non perforated




Boden gelocht
Bottom perforated



	außen · outside	innen · inside	
99-105-62	600 x 400 x 120 mm	550 x 350 x 100 mm	2189 g
99-105-68	600 x 400 x 180 mm	550 x 350 x 160 mm	2375 g

	außen · outside	innen · inside	
99-106-62	600 x 400 x 120 mm	550 x 350 x 100 mm	2431 g
99-106-68	600 x 400 x 180 mm	550 x 350 x 160 mm	2911 g

Deckel gelocht
Lid perforated



600 x 400 mm


Grey	99-100-51
Yellow	99-100-52
Red	99-100-53
Blue	99-100-54
Green	99-100-55
Black	99-100-56

Gewicht · weight
1668 g

Andere Container Other Containers



Boden ungelocht
Bottom non perforated



	99-310-050
	99-320-050
	99-330-050
	99-340-050
	99-350-050
	99-360-050

Gewicht · weight
1285 g



außen · outside	285 x 280 x 55 mm
innen · inside	270 x 270 x 33 mm

Boden gelocht
Bottom perforated



	99-310-080
	99-320-080
	99-330-080
	99-340-080
	99-350-080
	99-360-080

Gewicht · weight
1377 g



außen · outside	285 x 280 x 85 mm
innen · inside	270 x 270 x 64 mm

	99-311-080
	99-321-080
	99-331-080
	99-341-080
	99-351-080
	99-361-080


Gewicht · weight
1529 g

Boden ungelocht
Bottom non perforated




	außen · outside	innen · inside	
99-305-050	285 x 280 x 55 mm	270 x 270 x 33 mm	611 g
99-305-080	285 x 280 x 85 mm	270 x 270 x 64 mm	703 g

Boden gelocht
Bottom perforated



	außen · outside	innen · inside	
99-306-080	285 x 280 x 85 mm	270 x 270 x 64 mm	857 g

Deckel gelocht
Lid perforated



	99-304-010
	99-304-020
	99-304-030
	99-304-040
	99-304-050
	99-304-060


Gewicht · weight
674 g




Spannrahmen · Drape Retainer
99-395-000 255 x 255 mm



Drahtsiebkorb · Wire Basket
99-390-030 255 x 245 x 30 mm
99-390-050 255 x 245 x 50 mm



Zwischenboden · Wire Base
99-393-000 245 x 255 x 30 mm 351 g




Silikonmatte · Silicone Mat
99-851-000 245 x 245 mm



Einschlagtuch · Drape
99-095-040 90 x 60 cm
99-095-050 90 x 80 cm

Boden ungelocht
Bottom non perforated



	99-410-040
	99-420-040
	99-430-040
	99-440-040
	99-450-040
	99-460-040

Gewicht · weight
510,5 g



außen · outside	300 x 140 x 40 mm
innen · inside	280 x 130 x 27 mm

Mini Container

Kleinssetcontainer · Mini Container

Boden gelocht
Bottom perforated



	99-411-040
	99-421-040
	99-431-040
	99-441-040
	99-451-040
	99-461-040

Gewicht · weight
589 g

	99-411-070
	99-421-070
	99-431-070
	99-441-070
	99-451-070
	99-461-070

Gewicht · weight
630 g



außen · outside	300 x 140 x 70 mm
innen · inside	280 x 130 x 55 mm

	99-410-070
	99-420-070
	99-430-070
	99-440-070
	99-450-070
	99-460-070

Gewicht · weight
567 g

	99-411-100
	99-421-100
	99-431-100
	99-441-100
	99-451-100
	99-461-100


Gewicht · weight
625 g



außen · outside	300 x 140 x 100 mm
innen · inside	280 x 130 x 84 mm

	99-410-100
	99-420-100
	99-430-100
	99-440-100
	99-450-100
	99-460-100

Gewicht · weight
625,5 g



Mini Siebkorb · Mini Wire basket		
99-490-050	265 x 120 x 45 mm	251 g




Boden ungelocht
Bottom non perforated




	außen · outside	innen · inside	
99- 405-040	300 x 140 x 40 mm	280 x 130 x 27 mm	190 g
99-405-070	300 x 140 x 70 mm	280 x 130 x 55 mm	241 g
99-405-100	300 x 140 x 100 mm	280 x 130 x 84 mm	305 g



Boden gelocht
Bottom perforated



	außen · outside	innen · inside	
99-406-040	300 x 140 x 40 mm	280 x 130 x 27 mm	258 g
99-406-070	300 x 140 x 70 mm	280 x 130 x 55 mm	304 g
99-406-100	300 x 140 x 100 mm	280 x 130 x 84 mm	326 g



Deckel gelocht
Lid perforated

	99-400-010
	99-400-020
	99-400-030
	99-400-040
	99-400-050
	99-400-060

Gewicht · weight
325,5 g




Spannrahmen · Drape Retainer	
99-495-000	280 x 125 mm



Silikonmatte · Silicone Mat	
99-850-000	270 x 125 mm



Einschlagtuch · Drape	
99-095-010	60 x 40 cm



99-510-040

Gewicht · weight
783 g



outside	310 x 190 x 40 mm
inside	280 x 182 x 35 mm

99-510-060

Gewicht · weight
862 g




outside	310 x 190 x 65 mm
inside	280 x 182 x 55 mm

99-510-130



Gewicht · weight
1022 g



outside	310 x 190 x 130 mm
inside	280 x 182 x 124 mm



	outside	inside	
99-506-040	310 x 190 x 40 mm	280 x 182 x 35 mm	341 g
99-506-060	310 x 190 x 65 mm	280 x 182 x 55 mm	415 g
99-506-130	310 x 190 x 130 mm	280 x 182 x 124 mm	570 g



Deckel · Lid		
99-500-010	310 x 190 mm	417 g
99-500-020	310 x 190 mm	418 g



Spannrahmen · Drape Retainer	
99-595-000	310 x 190mm



Einschlagtuch · Drape	
99-095-010	60 x 40 cm
99-095-020	65 x 65 cm
99-095-030	85 x 85 cm



	99-612-070
	99-622-070
	99-632-070
	99-642-070
	99-652-070
	99-662-070

außen • outside	500 x 155 x 75 mm
innen • inside	495 x 155 x 75 mm

Gewicht • weight
654 g

Deckel gelocht • Boden gelocht (kein Filter)
Lid perforated • Bottom perforated (non Filter)




außen • outside	300 x 138 x 65 mm
innen • inside	295 x 138 x 65 mm

Gewicht • weight
368 g

	99-610-060
	99-620-060
	99-630-060
	99-640-060
	99-650-060
	99-660-060

Deckel gelocht · Boden gelocht (kein Filter)
Lid perforated · Bottom perforated (non Filter)



	99-600-020
--	------------

outside	450 x 70 x 70 mm
inside	450 x 70 x 60 mm

Gewicht · weight
302 g




	99-600-010
--	------------

outside	250 x 60 x 50 mm
inside	250 x 55 x 50 mm

Gewicht · weight
152 g



Deckel und Boden ungelocht
Lid and Bottom non perforated



	99-315-100
	99-325-100
	99-335-100
	99-345-100
	99-355-100
	99-365-100

außen · outside	285 x 280 x 100 mm
innen · inside	262 x 260 x 82 mm

Gewicht · weight
1677 g



	99-315-150
	99-325-150
	99-335-150
	99-345-150
	99-355-150
	99-365-150

außen · outside	285 x 280 x 150 mm
innen · inside	262 x 260 x 132 mm

Gewicht · weight
1929 g



	99-115-150
	99-125-150
	99-135-150
	99-145-150
	99-155-150
	99-165-150

außen · outside	580 x 280 x 150 mm
innen · inside	548 x 258 x 130 mm

Gewicht · weight
2883 g



Boden ungelocht
Bottom non perforated



99-105-150	580 x 280 x 150 mm	2108 g
------------	--------------------	--------



99-305-100	285 x 280 x 150 mm	1096 g
99-305-150	285 x 280 x 150 mm	1327 g

Deckel ungelocht
Lid non perforated



285 x 280 mm

	99-302-010
	99-302-020
	99-302-030
	99-302-040
	99-302-050
	99-302-060

Gewicht · weight

458 g

Deckel ungelocht
Lid non perforated



580 x 280 mm

	99-102-010
	99-102-020
	99-102-030
	99-102-040
	99-102-050
	99-102-060

Gewicht · weight

875 g



	99-882-000
--	------------

Brush container- For 20 brushes
240 x 275 x 70 mm
Stainless steel
Brushes are taken out by rotating the handle under the container-
Containers are supplied together with wall mounting fixtures

	99-880-100
	99-880-110
	99-880-120
	99-880-130
	99-880-140

Brush container- For 20 brushes
240 x 275 x 70 mm
Aluminium
Brushes are taken out by rotating the handle under the container-
Containers are supplied together with wall mounting fixtures




99-800-000
Disposable Paper Filter
VE 100 St-
Für Container
1/1, 3/4 und 1/2
Pack of 100
For 1/1, 3/4 and
Half Size Containers



99-800-010
Textilfilter
Für Container
1/1, 3/4 und 1/2
Single textile
For 1/1, 3/4 and
Half Size Containers



99-802-000
Disposable Paper Filter
Pack of 100
For Mini Containers
99-804-000
100 Pcs/Package
for Dental Containers



99-802-010
Textilfilter / Single textil
For Mini Containers
99-804-010
Textilfilter / Single textil
for Dental Containers




99-810-000
Pack of 100
Indikatorschilder · Paper label with indicator
For 1/1, 3/4 and Half Size Containers



99-812-000
Pack of 100
Indikatorschilder · Paper label with indicator
For Mini and Dental Containers



Filterhalterung · filter tray
99-809-000 For Paper Filter



99-801-000
Spender für Einmalfilter
Filter Dispenser



99-819-000
Sicherheitssiegel
Plastic Security Seal
Pack of 100




99-829-000
Sortiermagazin für Kennzeichnungsschilder
Dispenser for identification labels




	99-820-010
Spinal	99-820-020
E.N.T	99-820-030
General	99-820-040
Neurosurgery	99-820-050
Micro	99-820-060
	99-820-070
	99-820-080

	99-820-010B
	99-820-020B
	99-820-030B
	99-820-040B
	99-820-050B
	99-820-060B
	99-820-070B
	99-820-080B

Kennzeichnungsschilder mit oder ohne
Beschriftung, max- 16 Buchstaben
Coloured identification labels- Printed to your specificati-
on, max-16 letters





270 x 125 mm



245 x 245 mm

Size	gelocht grün	gelocht blau	ungelocht grün	ungelocht blau
270 x 125 mm	99-850-001G	99-850-001B	99-850-000G	99-850-000B
245 x 245 mm	99-851-001G	99-851-001B	99-851-000G	99-851-000B
220 x 230 mm	99-852-001G	99-852-001B	99-852-000G	99-852-000B
380 x 230 mm	99-854-001G	99-854-001B	99-854-000G	99-854-000B
440 x 230 mm	99-854-201G	99-854-201B	99-854-200G	99-854-200B
520 x 230 mm	99-854-301G	99-854-301B	99-854-300G	99-854-300B



99-849-200
24 x 28 x 25
ø 7 - 12 mm



99-849-210
24 x 28 x 29
ø 12 - 16 mm



99-849-220
24 x 28 x 31
ø 16 - 26 mm



99-849-230
24 x 28 x 40
ø 26 - 36 mm





99-849-240
24 x 28 x 44
ø 36 - 45 mm



99-849-000
50 x 20 x 12 mm



99-849-010
65 x 20 x 20 mm






99-830-010
Haltestift für Trennbleche
Holding Pin for Partition Sheets



99-830-000
Haltestifte für Instrumente
Holding Pin for Wire Baskets



Trennbleche · Sheets	
99-831-050	50 x 20 mm
99-831-130	130 x 20 mm
99-831-220	225 x 20 mm
99-831-460	460 x 20 mm



Schiene für Gleitschieber
Sliding Rails with Securing Screws
99--841-000 230 mm
99--841-010 470 mm



99-842-000 230 mm
99-842-010 470 mm
Anschlagrahmen
Limiting Bar with Securing Screws




99-840-000
Gleitschieber
Divider, to be used
with 99-841-000



99-845-000
Spiralfeder Haltevorrichtung
Spiral Holding




99-847-000
Doppel-Spiralfeder Haltevorrichtung
Double Spiral Holding



99-800-020
Verwendungsdauer 1200 Sterilisationszyklen
für 1/1, 3/4 und 1/2 Container
Life for 1200 Sterilization cycles
for 1/1, 3/4 and 1/2 Container



99-804-020
Verwendungsdauer 1200 Sterilisationszyklen
für Dental Container
Life for 1200 Sterilization cycles
for Dental Container



99-802-020
Verwendungsdauer 1200 Sterilisationszyklen
für Mini Container
Life for 1200 Sterilization cycles
for Mini Container

Siebkörbe




Anhand dieser Übersicht sehen Sie die verwendeten Standard Siebkorbgrößen je Container Modell.

This overview shows the standard used basket sizes for each container model.


Container Modell		Standard Siebgrößen
1/1 Container	580 x 280 mm	540 x 250 mm
		480 x 250 mm
3/4 Container	465 x 280 mm	405 x 250 mm
1/2 Container	285 x 280 mm	250 x 240 mm
Mini Container	300 X 140 mm	240 x 120 mm
		250 x 120 mm
		265 x 120 mm
Dental Container	310 x 190 mm	267 x 170 mm
		280 x 170 mm

Einschlagtücher / Drapes

Item Code	Artbez.
99-095-000	OP-Tuch, Gewebebreite 160cm, grün
99-095-010	Tuch 60x40 cm, grün
99-095-020	Tuch 60x60 cm, grün
99-095-030	Tuch 85x85 cm, grün
99-095-040	Tuch 90x60 cm, grün
99-095-050	Tuch 90x80 cm, grün
99-095-060	Tuch 130x45 cm, grün
99-095-070	Tuch 130x60 cm, grün
99-095-080	Tuch 140x100 cm, grün
99-095-090	Tuch 185x45 cm, grün
99-095-100	Tuch 185x60 cm, grün



99-095-01B	Tuch 60x40 cm, blau
99-095-02B	Tuch 60x60 cm, blau
99-095-03B	Tuch 85x85 cm, blau
99-095-04B	Tuch 90x60 cm, blau
99-095-05B	Tuch 90x80 cm, blau
99-095-06B	Tuch 130x45 cm, blau
99-095-07B	Tuch 130x60 cm, blau
99-095-08B	Tuch 140x100 cm, blau
99-095-09B	Tuch 185x45 cm, blau
99-095-10B	Tuch 185x60 cm, blau




Siebkörbe aus Edelstahl mit oder ohne Deckel

Sondergrößen auf Anfrage


Baskets from stainless Steel with or without lid

Special sizes on inquiry



Item Code	Size in mm	Cont.
98-590-050	Siebkorb aus Edelstahlblech 267 x 173 x 34 mm	Dental
98-490-050	Siebkorb aus Edelstahlblech 271 x 117 x 50 mm	Mini
98-390-030	Siebkorb aus Edelstahlblech 243 x 253 x 30 mm	1/2
98-390-050	Siebkorb aus Edelstahlblech 243 x 253 x 50 mm	1/2
98-390-070	Siebkorb aus Edelstahlblech 243 x 253 x 70 mm	1/2
98-390-100	Siebkorb aus Edelstahlblech 243 x 253 x 100 mm	1/2
98-290-030	Siebkorb aus Edelstahlblech 405 x 253 x 30 mm	3/4
98-290-050	Siebkorb aus Edelstahlblech 405 x 253 x 50 mm	3/4
98-290-070	Siebkorb aus Edelstahlblech 405 x 253 x 70 mm	3/4
98-290-100	Siebkorb aus Edelstahlblech 405 x 253 x 100 mm	3/4
98-190-030	Siebkorb aus Edelstahlblech 485 x 253 x 30 mm	1/1
98-190-050	Siebkorb aus Edelstahlblech 485 x 253 x 50 mm	1/1
98-190-070	Siebkorb aus Edelstahlblech 485 x 253 x 70 mm	1/1
98-190-100	Siebkorb aus Edelstahlblech 485 x 253 x 100 mm	1/1
98-191-030	Siebkorb aus Edelstahlblech 540 x 253 x 30 mm	1/1
98-191-050	Siebkorb aus Edelstahlblech 540 x 253 x 50 mm	1/1
98-191-070	Siebkorb aus Edelstahlblech 540 x 253 x 70 mm	1/1
98-191-100	Siebkorb aus Edelstahlblech 540 x 253 x 100 mm	1/1

98-590-000	Deckel zu Siebkorb 265 x 175 mm mit Verschuß	Dental
98-490-000	Deckel zu Siebkorb 271 x 117 mm mit Verschuß	Mini
98-390-000	Deckel zu Siebkorb 243 x 253 mm mit Verschuß	1/2
98-290-000	Deckel zu Siebkorb 405 x 253 mm mit Verschuß	3/4
98-190-000	Deckel zu Siebkorb 485 x 253 mm mit Verschuß	1/1
98-191-000	Deckel zu Siebkorb 540 x 253 mm mit Verschuß	1/1



Zwischenboden für Edelstahlsiebkorb, für alle Siebkorbgrößen lieferbar

Wire base for baskets, available for all sizes

Universal Silicon Holders



Item Code	Sizes in mm
99-8950	250x35x10



Item Code	Sizes in mm
99-8951	250x35x10




Silikonlagerungen

Silikonlagerungen schützen Ihre Instrumente während der Sterilisation, dem Transport und der Lagerung. Die Einteilung Ihres Sets ist flexibel- Wir gestalten das Instrumentarium im Siebkorb nach Ihren Wünschen-

Silicone fixation

The silicone fixation protects your instruments during transportation, sterilisation and storage- The partitioning of the instruments is flexible and can be adapted to suit your specific requirements-



Our side perforated baskets are:

- Manufactured in stainless steel grade AISI-304 .
- Electro-polished.
- Burr free to protect, wraps and eliminate user injury.
- Smooth round Edges.
- Side wall perforated to protect sharp instruments edges and tips from damage. The perforation hole size is 2.0 mm.
- Meet all requirements of automatic washing and sterilization of all kinds of instruments.
- Separate lids are available with single, double and detention frame.
- Inner mounted drop handles.
- Top and bottom frame diameter 5.0 mm.
- Base available in wire mesh with opening square 6.0 mm, wire diameter 1.2 mm.
- Stacking frame can also be available on request.



Item Code	Sizes in mm	Cont.
99-8610	240X120x30	Mini
99-8611	240X120x50	Mini
99-8612	240X120x60	Mini
99-8608	240X120x70	Mini
99-8609	240X120x100	Mini

99-8607	240X250x30	1/2
99-8613	240X250x50	1/2
99-8615	240X250x60	1/2
99-8616	240X250x70	1/2
99-8617	240X250x100	1/2

99-8656	405X250x30	3/4
99-8657	405X250x50	3/4
99-8658	405X250x60	3/4
99-8659	405X250x70	3/4
99-8660	405X250x100	3/4

99-8606	480X250x30	1/1
99-8618	480X250x50	1/1
99-8619	480X250x60	1/1
99-8620	480X250x70	1/1
99-8621	480X250x100	1/1

Item Code	Sizes in mm	Cont.
99-8661	540x250x30	1/1
99-8662	540x250x50	1/1
99-8663	540x250x60	1/1
99-8664	540x250x70	1/1
99-8665	540x250x100	1/1

99-8622	300x300x50
99-8623	300x300x60


99-8624	470x400x60
---------	------------

99-8625	480x480x50
99-8626	480x480x60


99-8628	500x300x50
---------	------------

99-8632	550x430x50
---------	------------

99-8634	600x300x60
---------	------------



Smooth Round Edges



Inner Mounted Drop Handles

Many more sizes available on request

Lid for perforated baskets (Single Frame)


Item Code	Sizes in mm	Cont.
99-8635	240X120	Mini
99-8636	240X250	1/2
99-8642	405x250	3/4
99-8637	480X250	1/1
99-8643	540x250	1/1
99-8638	300X300	
99-8639	480X480	
99-8640	500X300	
99-8641	600X300	



Many more sizes available on request

Lid for perforated baskets (Double Frame)


Item Code	Sizes in mm	Cont.
99-8635D	240x120	Mini
99-8636D	240X250	1/2
99-8642D	405x250	3/4
99-8637D	480X250	1/1
99-8643D	540x250	1/1
99-8638D	300X300	
99-8639D	480X480	
99-8640D	500X300	
99-8641D	600X300	




Many more sizes available on request

Lid for perforated baskets (Detention Frame)


Item Code	Sizes in mm	Cont.
99-8635DT	240X120	Mini
99-8636DT	240X250	1/2
99-8642DT	405x250	3/4
99-8637DT	480X250	1/1
99-8643DT	540x250	1/1
99-8638DT	300X300	
99-8639DT	480X480	
99-8640DT	500X300	
99-8641DT	600X300	



Many more sizes available on request



Flat base (Inner Side)




Inner mounted drop handle



Flat base (Outer Side)



Smooth round edges



Item Code	Sizes in mm	Cont.
99-8545	120x250x50	Mini
99-8546	120x250x60	Mini
99-8547	120x250x70	Mini
99-8548	120x250x100	Mini
99-8549	240x250x50	1/2
99-8550	240x250x60	1/2
99-8551	240x250x70	1/2
99-8552	240x250x100	1/2
99-8553	405x250x50	3/4
99-8554	405x250x60	3/4
99-8555	405x250x70	3/4
99-8556	405x250x100	3/4
99-8557	480x250x50	1/1
99-8558	480x250x60	1/1
99-8559	480x250x70	1/1
99-8560	480x250x100	1/1

Item Code	Sizes in mm	Cont.
99-8561	540x250x50	1/1
99-8562	540x250x60	1/1
99-8563	540x250x70	1/1
99-8564	540x250x100	1/1
99-8565	480x360x50	
99-8566	480x360x60	
99-8567	480x360x70	
99-8568	480x360x100	
99-8569	480x480x50	
99-8570	480x480x60	
99-8571	480x480x70	
99-8572	480x480x100	



Universal Silicon Holders



Item Code	Sizes in mm
99-8950	250x35x10




Item Code	Sizes in mm
99-8951	250x35x10

Side Perforated Basket Flat Base

Our side perforated baskets are:

- Manufactured in stainless steel grade AISI-304 .
- Electro-polished.
- Burr free to protect, wraps and eliminate user injury.
- Smooth round Edges.
- Side wall perforated to protect sharp instruments
- edges and tips from damage.
- The perforation hole size is 2.0 mm.
- Meet all requirements of automatic washing and sterilization
- of all kinds of instruments.
- Separate lids are available with double frame.
- Inner mounted drop handles.
- Top frame diameter 5.0 mm.
- Flat base basket available in welded mesh with opening
- square 6.0mm, wire diameter 1.2 mm.
- Stacking frame can also be available on request.



Lid for perforated baskets (Double Frame)


Item Code	Sizes in mm	Cont.
99-8582	120x250	Mini
99-8583	240x250	1/2
99-8584	405x250	3/4
99-8585	480x250	1/1
99-8586	540x250	1/1
99-8587	480x360	
99-8588	480x480	

Many more sizes available on request




Standard Wire Mesh Baskets

- Our wire mesh baskets are:
- Manufactured in stainless steel grade AISI-304.
- Electro-polished.
- Burr free to protect, wraps and eliminate user injury.
- Meet all requirements of automatic washing and sterilization of all kinds of instruments.
- Separate lids are available with single, double and detention frame.
- Available in drop and folding handles.
- Top and bottom frame diameter 5.0 mm.
- Wire mesh with opening square 4.5 mm,
- wire diameter 1.00 mm.
- Baskets are stackable.




Inner Mounted Drop Handles



Corner Edge



Bottom Stacking Frame



Stackability

Item Code	Sizes in mm	Cont.
99-8401	105x240x30	
99-8402	105x240x50	
99-8403	105x240x60	
99-8404	105x240x70	
99-8405	105x240x100	
99-8406	240x240x30	
99-8407	240x240x50	
99-8408	240x240x60	
99-8409	240x240x70	
99-8410	240x240x100	
99-8411	240x250x30	1/2
99-8412	240x250x50	1/2
99-8413	240x250x60	1/2
99-8414	240x250x70	1/2
99-8415	240x250x100	1/2
99-8416	250x250x30	
99-8417	250x250x50	
99-8418	250x250x60	
99-8419	250x250x70	
99-8420	250x250x100	

Item Code	Sizes in mm	Cont.
99-8421	405X250X30	3/4
99-8422	405X250X50	3/4
99-8423	405X250X60	3/4
99-8424	405X250X70	3/4
99-8425	405X250X100	3/4
99-8426	480X240X30	
99-8427	480X240X50	
99-8428	480X240X60	
99-8429	480X240X70	
99-8430	480X240X100	
99-8431	480X250X30	1/1
99-8432	480X250X50	1/1
99-8433	480X250X60	1/1
99-8434	480X250X70	1/1
99-8435	480X250X100	1/1
99-8436	480X480X30	
99-8437	480X480X50	
99-8438	480X480X60	
99-8439	480X480X70	
99-8440	480X480X100	

Item Code	Sizes in mm	Cont.
99-8441	540X250X30	1/1
99-8442	540X250X50	1/1
99-8443	540X250X60	1/1
99-8444	540X250X70	1/1
99-8445	540X250X100	1/1

Many more sizes available on request

Universal Silicon Holders



Item Code	Sizes in mm
99-8950	250x35x10



Item Code	Sizes in mm
99-8951	250x35x10

Lids for Wire Mesh Baskets (Single Frame)

Item Code	Sizes in mm	Cont.
99-8455	105X240	
99-8456	240X240	
99-8457	240X250	1/2
99-8458	250X250	
99-8459	405X250	3/4
99-8460	480X240	
99-8461	480X250	1/1
99-8462	480X480	
99-8463	540X250	1/1


Many more sizes available on request



Lids for Wire Mesh Baskets (Double Frame)

Item Code	Sizes in mm	Cont.
99-8455D	105X240	
99-8456D	240X240	
99-8457D	240X250	1/2
99-8458D	250X250	
99-8459D	405X250	3/4
99-8460D	480X240	
99-8461D	480X250	1/1
99-8462D	480X480	
99-8463D	540X250	1/1


Many more sizes available on request



Lids for Wire Mesh Baskets (Detention Frame)

Item Code	Sizes in mm	Cont.
99-8455DT	105X240	
99-8456DT	240X240	
99-8457DT	240X250	1/2
99-8458DT	250X250	
99-8459DT	405X250	3/4
99-8460DT	480X240	
99-8461DT	480X250	1/1
99-8462DT	480X480	
99-8463DT	540X250	1/1

Many more sizes available on request






Crimped Wire Mesh Baskets

Our crimped wire mesh baskets are:

- Manufactured in stainless steel grade AISI-304.
- Electro-polished.
- Burr free to protect, wraps and eliminate user injury.
- Meet all requirements of automatic washing and sterilization of all kinds of instruments.
- Separate lids are available with single, double and detention frame.
- Available in tilted and drop handles.
- Top and bottom frame diameter 5.0 mm.
- Crimped mesh with opening square 5.0 mm,
- wire diameter 0.9 mm.
- Baskets are stackable.



Stackability

Bottom Stacking Frame

Tilted Handle

Item Code	Sizes in mm	Cont.	Item Code	Sizes in mm	Cont.	Item Code	Sizes in mm	Cont.
99-8401CW	105X240X30		99-8421CW	405X250X30	3/4	99-8441CW	540X250X30	1/1
99-8402CW	105X240X50		99-8422CW	405X250X50	3/4	99-8442CW	540X250X50	1/1
99-8403CW	105X240X60		99-8423CW	405X250X60	3/4	99-8443CW	540X250X60	1/1
99-8404CW	105X240X70		99-8424CW	405X250X70	3/4	99-8444CW	540X250X70	1/1
99-8405CW	105X240X100		99-8425CW	405X250X100	3/4	99-8445CW	540X250X100	1/1
99-8406CW	240X240X30		99-8426CW	480X240X30				
99-8407CW	240X240X50		99-8427CW	480X240X50				
99-8408CW	240X240X60		99-8428CW	480X240X60				
99-8409CW	240X240X70		99-8429CW	480X240X70				
99-8410CW	240X240X100		99-8430CW	480X240X100				
99-8411CW	240X250X30	1/2	99-8431CW	480X250X30	1/1			
99-8412CW	240X250X50	1/2	99-8432CW	480X250X50	1/1			
99-8413CW	240X250X60	1/2	99-8433CW	480X250X60	1/1			
99-8414CW	240X250X70	1/2	99-8434CW	480X250X70	1/1			
99-8415CW	240X250X100	1/2	99-8435CW	480X250X100	1/1			
99-8416CW	250X250X30		99-8436CW	480X480X30				
99-8417CW	250X250X50		99-8437CW	480X480X50				
99-8418CW	250X250X60		99-8438CW	480X480X60				
99-8419CW	250X250X70		99-8439CW	480X480X70				
99-8420CW	250X250X100		99-8440CW	480X480X100				

Many more sizes
available on request

Universal Silicon Holders



Item Code	Sizes in mm
99-8950	250x35x10



Item Code	Sizes in mm
99-8951	250x35x10

Lid for Crimped Wire Mesh Baskets (Single Frame)

Item Code	Sizes in mm	Cont.
99-8455CW	105X240	
99-8456CW	240X240	
99-8457CW	240X250	1/2
99-8458CW	250X250	
99-8459CW	405X250	3/4
99-8460CW	480X240	
99-8461CW	480X250	1/1
99-8462CW	480X480	
99-8463CW	540X250	1/1


Many more sizes available on request



Lid for Crimped Wire Mesh Baskets (Double Frame)

Item Code	Sizes in mm	Cont.
99-8455DCW	105X240	
99-8456DCW	240X240	
99-8457DCW	240X250	1/2
99-8458DCW	250X250	
99-8459DCW	405X250	3/4
99-8460DCW	480X240	
99-8461DCW	480X250	1/1
99-8462DCW	480X480	
99-8463DCW	540X250	1/1


Many more sizes available on request



Lid for Crimped Wire Mesh Baskets (Detention Frame)

Item Code	Sizes in mm	Cont.
99-8455DTCW	105X240	
99-8456DTCW	240X240	
99-8457DTCW	240X250	1/2
99-8458DTCW	250X250	
99-8459DTCW	405X250	3/4
99-8460DTCW	480X240	
99-8461DTCW	480X250	1/1
99-8462DTCW	480X480	
99-8463DTCW	540X250	1/1

Many more sizes available on request







Classic Crimped Wire Mesh Baskets

Our Classic crimped wire mesh baskets are:

- Manufactured in stainless steel grade AISI-304.
- Electro-polished.
- Burr free to protect, wraps and eliminate user injury.
- Meet all requirements of automatic washing and sterilization
- of all kinds of instruments.
- Straight and smooth round corners.
- Available in tilted and drop handles.
- Top and bottom frame diameter 5.0 mm.
- Crimped mesh with opening square 6.0 mm,
- wire diameter 1.0 mm.



Smooth and Straight Corner



Tilted Handle

Item Code	Sizes in mm	Cont.	Item Code	Sizes in mm	Cont.	Item Code	Sizes in mm	Cont.								
99-8401ER	105X240X30		99-8421ER	405X250X30	3/4	99-8441ER	540X250X30	1/1								
99-8402ER	105X240X50		99-8422ER	405X250X50	3/4	99-8442ER	540X250X50	1/1								
99-8403ER	105X240X60		99-8423ER	405X250X60	3/4	99-8443ER	540X250X60	1/1								
99-8404ER	105X240X70		99-8424ER	405X250X70	3/4	99-8444ER	540X250X70	1/1								
99-8405ER	105X240X100		99-8425ER	405X250X100	3/4	99-8445ER	540X250X100	1/1								
99-8406ER	240X240X30		99-8426ER	480X240X30	Many more sizes available on request	Universal Silicon Holders  <table><tr><th>Item Code</th><th>Sizes in mm</th></tr><tr><td>99-8950</td><td>250x35x10</td></tr></table>  <table><tr><th>Item Code</th><th>Sizes in mm</th></tr><tr><td>99-8951</td><td>250x35x10</td></tr></table>			Item Code	Sizes in mm	99-8950	250x35x10	Item Code	Sizes in mm	99-8951	250x35x10
Item Code	Sizes in mm															
99-8950	250x35x10															
Item Code	Sizes in mm															
99-8951	250x35x10															
99-8407ER	240X240X50	99-8427ER	480X240X50													
99-8408ER	240X240X60	99-8428ER	480X240X60													
99-8409ER	240X240X70	99-8429ER	480X240X70													
99-8410ER	240X240X100	99-8430ER	480X240X100													
99-8411ER	240X250X30	99-8431ER	480X250X30	1/1												
99-8412ER	240X250X50	99-8432ER	480X250X50	1/1												
99-8413ER	240X250X60	99-8433ER	480X250X60	1/1												
99-8414ER	240X250X70	99-8434ER	480X250X70	1/1												
99-8415ER	240X250X100	99-8435ER	480X250X100	1/1												
99-8416ER	250X250X30	99-8436ER	480X480X30													
99-8417ER	250X250X50	99-8437ER	480X480X50													
99-8418ER	250X250X60	99-8438ER	480X480X60													
99-8419ER	250X250X70	99-8439ER	480X480X70													
99-8420ER	250X250X100	99-8440ER	480X480X100													



Fine Mesh Baskets With Lid


Our Fine Mesh Baskets With Lid Are:

- Manufactured in stainless steel grade AISI-304.
- Electro-polished.
- Burr free to protect, wraps and eliminate user injury. Smooth round Edges.
- Top and lid frame wire diameter 3.5 mm and 04.0 mm (according to the size of baskets).
- Fine mesh with opening square 1.6 mm and wire diameter 0.7 mm
- Wire Locking Clip.

Item Code	Sizes in mm
99-8343	40x40x20
99-8344	70x70x30
99-8345	80x40x20
99-8346	80x40x30
99-8347	80x80x34
99-8348	80x80x40
99-8349	90x90x50
99-8350	150x110x20
99-8351	160x40x20
99-8352	105x70x25
99-8353	141x50x40
99-8354	150x50x35
99-8355	180x130x20
99-8356	190x85x100
99-8357	200x100x35
99-8358	220x140x50
99-8359	210x170x20
99-8360	235x155x40

Item Code	Sizes in mm	Cont.
99-8361	240x90x40	
99-8362	240x120x35	Mini
99-8363	275x175x35	
99-8364	280x170x60	Dental
99-8365	250x200x20	
99-8366	300x180x20	
99-8367	200x110x40	
99-8368	225x110x50	
99-8369	235x120x60	
99-8370	270x135x40	
99-8371	280x145x50	
99-8372	400x250x50	


Many more sizes
available on request



Micro Fine Mesh Baskets

Our Micro Fine Mesh Baskets With Lid Are:

- Manufacture in stainless steel grade AISI-304.
- Top and lid frame wire diameter 3.5 mm
- Micro fine mesh with opening square 1.0 mm
- and wire diameter 0.5 mm
- Wire locking Clip.
- Electrolyte Polish.



Item Code	Sizes in mm
99-8343MM	40X40x20
99-8344MM	70x70x30
99-8345MM	80x40x20
99-8346MM	80x40x30
99-8347MM	80x80x34
99-8348MM	80x80x40
99-8349MM	90x90x50
99-8350MM	150x110x20
99-8351MM	160x40x20
99-8352MM	105x70x25
99-8353MM	141x50x40
99-8354MM	150x50x35
99-8355MM	180x130x20

Many more sizes available on request

Fine Mesh Baskets with Unique Twist-Off Lid

Our Fine Mesh Baskets With Lid Are:

- Manufacture in High Quality stainless steel grade AISI-304.
- Top and lid frame wire diameter 3.5 mm and 4.0 mm (According to the size of Basket)
- Fine mesh with opening square 1.6 mm
- and wire diameter 0.7 mm
- Burr-Free and Electrolyte Polish.
- Innovative twist-off Lids are secure and easy to use



Item Code	Sizes in mm
99-8344RL	70x70x30
99-8347RL	80x80x34
99-8348RL	80x80x40
99-8349RL	90x90x50
99-8350RL	150x110x20
99-8352RL	105x70x25
99-8353RL	141x50x40
99-8354RL	150x50x35
99-8355RL	180x130x20
99-8356RL	190x85x100
99-8357RL	200x100x35
99-8358RL	220x140x50
99-8359RL	210x170x20
99-8360RL	235x155x40

Item Code	Sizes in mm	Cont.
99-8361RL	240x90x40	
99-8362RL	240x120x35	Mini
99-8363RL	275x175x35	
99-8364RL	280x170x60	Dental
99-8365RL	250x200x20	
99-8366RL	300x180x20	
99-8367RL	200x110x40	
99-8368RL	225x110x50	
99-8369RL	235x120x60	
99-8370RL	270x135x40	
99-8371RL	280x145x50	
99-8372RL	400x250x50	

Many more sizes available on request

Micro Fine Mesh Baskets with Unique Twist-Off Lid

Our Micro Fine Mesh Baskets With Lid Are:

- Manufacture in High Quality stainless steel grade AISI-304.
- Top and lid frame wire diameter 3.5 mm
- Micro Fine mesh with opening square 1.0 mm
- and wire diameter 0.5 mm
- Burr-Free and Electrolyte Polish.
- Innovative twist-off Lids are secure and easy to use




Item Code	Sizes in mm
99-8344MMRL	70x70x30
99-8347MMRL	80x80x34
99-8348MMRL	80x80x40
99-8349MMRL	90x90x50
99-8350MMRL	150x110x20
99-8352MMRL	105x70x25
99-8353MMRL	141x50x40
99-8354MMRL	150x50x35
99-8355MMRL	180x130x20

Many more sizes available on request

Our Full perforated trays are:


- Manufacture in stainless steel grade AISI-304.
- Electro Polished.
- Top and lid frame rod diameter 3.5 mm and 4.0 mm.
- Full wall perforation to protect sharp instruments
- edges and tips from damage. The perforation hole size is 2.0mm
- Smooth round edges.
- Sheet locking clip.
- Hinged Lids.



Item Code	Sizes in mm
99-6343	40X40X20mm
99-6344	70x70x30mm
99-6345	80x40x20mm
99-6346	80x40x30mm
99-6347	80x80x34mm
99-6348	80x80x40mm
99-6349	90x90x50mm
99-6350	150x110x20mm
99-6351	160x40x20mm
99-6352	105x40x20mm
99-6353	141x50x40mm
99-6354	150x50x35mm
99-6355	190x85x35mm
99-6356	200x100x35mm
99-6357	240x90x40mm
99-6358	200x110x40mm
99-6359	225x110x50mm


Many more sizes available on request

- Bottom frame fabricated to make space between instruments and working surface.
- Separate top cover available.
- Stackable.
- Top and bottom frame wire diameter 3.5 mm
- Grid made with wire diameter 1.5 mm, opening square 6.0 mm
- Burr free smooth finish.
- Electrolyte Polish.



WIRE NET BASKET

Item Code	Sizes in mm
99-9000	250x170x70mm
99-9001	340x250x70mm
99-9002	450x340x70mm



LID FOR WIRE NET BASKET

Item Code	Sizes in mm
99-9003	250x170mm
99-9004	340x250mm
99-9005	450x340mm

- The standard baskets can be used for transportation, storage
- for all types of autoclaves. STE Sterile Goods Baskets means
- the baskets made with continuous wire which provide maximum
- protection to user from being injured. The baskets are stackable
- and can also be mounted on the cart or wall with hooks.





Item Code	Sizes in mm	
99-9010	575x280x135	Standard
99-9011	575x280x265	Standard
99-9012	575x280x265	Front access



- Upper frame wire diameter 6.0 mm
- Frame wire diameter 3.0 mm
- Loops wire diameter 2.5 mm
- Base frame wire diameter 5.0 mm
- Electrolyte Polish and burr free smooth finish.

- The SPRI baskets are both stackable and nestable.
- Stackable when the handles are folded into basket.
- Nestable when the handles are folded out.
- When the baskets are empty, hold minimum storage space in nesting position.
- Upper and lower frame wire diameter 5.0mm.
- Loops wire diameter 3.5mm
- Electrolyte polish and smooth finish.





**SPRI BASKETS
(Double Frame
with Folding Handle)**

Item Code	Sizes in mm
99-9110	585x395x100
99-9111	585x395x195




**SPRI BASKETS
(With Single Frame)**

Item Code	Sizes in mm
99-9030	600x300x40
99-9031	600x300x50
99-9032	600x300x100
99-9033	600x300x130
99-9034	600x300x150
99-9035	600x300x180
99-9036	600x300x200
99-9037	600x400x100
99-9038	600x400x130
99-9039	600x400x150
99-9040	600x400x180
99-9041	600x400x200




- This type of baskets are also both nestable and stackable.
- Stackable when placed the basket on each other with opposite position.
- Nestable when placed in same position.
- Upper and lower frame wire diameter 6.0 mm.
- Loops wire diameter 4.0 mm.
- Electrolyte Polish.

Item Code	Sizes in mm
99-9060	600x300
99-9061	600x400




Wire diameter 6.0mm and 2.5mm



Wire diameter 6.0mm and 3.0mm

Item Code	Sizes in mm
99-9062	600x400x44

Item Code	Sizes in mm
99-9063	600x400x44




Wire diameter 6.0mm and 3.0mm

Item Code	Sizes in mm	Scope
99-8046	640x150x77	1

Small Lens Baskets




Locking Clip



Arthroscope Basket




Item Code	Sizes in mm	Scope
99-8040	180X80X52	1
99-8041	290X80X52	1



Item Code	Sizes in mm	Scope
99-8042	290X80X52	2


Endoscope Baskets




Item Code	Sizes in mm	Scope
99-8044	460x80x52	2

Item Code	Sizes in mm	Scope
99-8043	460x80x52	1

Item Code	Sizes in mm	Scope
99-8045	670x80x52	1
99-8047	600x80x52	1
99-8048	700x80x52	1



Item Code	Sizes in mm	Scope
99-8049	290x160x52	2



Item Code	Sizes in mm	Scope
99-8050	290x160x52	2



Item Code	Sizes in mm	Scope
99-8051	460x160x52	2




Item Code	Sizes in mm	Scope
99-8052	460x160x52	2

- Top and Lid frame wire diameter 5.0mm.
- Loops wire diameter 2.0mm.
- Autoclaveable silicon inserts.
- Electrolyte polish and burr free finish.
- Locking Clips



Item Code	Sizes in mm	Scope
99-8053	460x160x52	2

Item Code	Sizes in mm	Scope
99-8054	310x150x77	2



Item Code	Sizes in mm	Scope
99-8066	460x225x52	2




Item Code	Sizes in mm	Scope
99-8067	460x225x52	2



Item Code	Sizes in mm	Scope
99-8068	580x160x52	2




Item Code	Sizes in mm	Scope
99-8069	580x160x52	2



Item Code	Sizes in mm	Cont..
99-822-000	240x250x38	1/2
99-822-001	480x250x38	1/1
99-822-002	400x250x38	



- With silicon inserts on bottom of baskets to accommodate
- micro instruments.
- Lids are provided with silicon strips to grip the instruments.
- Side walls are made of perforated sheet.
- User friendly locking handle.



Item Code	Sizes in mm
99-820-020	160X70X25
99-820-030	265X160X25
99-820-040	265X160X50






Spring Holding Dental Cassette

- Specially designed to minimize the contact area while locking of instruments inside the cassette.
- Each Instrument is held by separate spring.
- Further it can equally hold thicker and thinner instruments in one arrangement of cassette.

Item Code	Sizes in mm
99-8242	178X64x32
99-8243	178X114x32
99-8244	203X140x38
99-8245	203x280x38
99-8246	220x125x40
99-8247	250x220x40

Fine Mesh Dental Cassette

- Fine Mesh dental cassette is much workable when instruments are required
- to be cleaned by all around water flow.



Item Code	Sizes in mm
99-8148	240X90x40
99-8149	235x80x30
99-8150	235x155x40



Item Code	Sizes in mm
99-7010	180X60X22
99-7011	180X125X22



Item Code	Sizes in mm
99-7012	184x145x34
99-7013	184x74x34
99-7014	184x74x34





Item Code	Sizes in mm
99-7015	205x143x34



Item Code	Sizes in mm
99-7016	280x183x32

for 16 Instruments with partition plate



Item Code	Sizes in mm
99-7017	280x183x32

for 20 Instruments



Item Code	Sizes in mm
99-7018	280x205x32



Item Code	Sizes in mm
99-7019	280x183x32

Silicon inserts available in different colors



Green



orange



Yellow



Red



Blue



Micro cassettes are the best for sterilization of micro instruments and small clamps.

Item Code	Sizes in mm
99-7040	40x55x20
99-7041	74x45x20
99-7042	74x59x24



Clip-IN



Item Code	Sizes in mm
99-7051	185x82x18
99-7052	185x140x18

Fixing Clamp

Our Fixing Clamps are:

- Adjustable clamp length.
- Universally adaptable for all sizes of instruments.
- Easily dismantled for cleaning purpose.



Item Code	Sizes in mm
99-8141	60/220x220x39
99-8142	60x220x39

Item Code	Sizes in mm
99-7060	180x148x50



Item Code	Sizes in mm
99-7061	100x90x50

- Implantology cassettes are strongly recommended for the process of cleaning and sterilization of implant accessories.
- Different silicon plugs and holders are helpful in holding different burs and drills of 2 to 6.5 mm.
- Special silicon arrangement to hold ratchet.
- A separate removable compartments for various micro parts.
- Lids are also available without perforation.

Item Code	Sizes in mm
99-7062	100x90x50



Item Code	Sizes in mm
99-7063	150x85x80



Item Code	Sizes in mm
99-7064	150x85x80



Item Code	Sizes in mm
99-7065	80x55x40



Item Code	Sizes in mm
99-7066	80x55x40



Item Code	Sizes in mm
99-7067	80x55x40

Bur Clips



Item Code	Sizes in mm
99-7073	60x45x28 17 Holes

Item Code	Sizes in mm
99-7072	45x38x28 14 Holes

Item Code	Sizes in mm
99-7071	35x28x28 8 Holes

Bur Holder Stand



Item Code	Sizes in mm
99-7080	40x75x25 6 HOLES
99-7081	60x75x20 6 HOLES
99-7082	70x75x20 6 HOLES
99-7083	45x20x30 12 HOLES
99-7084	50x75x25 12 HOLES

Item Code	Sizes in mm
99-7090	45x20x30 4 HOLES
99-7091	35x25x30 6 HOLES
99-7092	45x25x30 8 HOLES
99-7093	63x25x30 12 HOLES
99-7094	85x25x30 16 HOLES
99-7095	110x35x34 30 HOLES

Instruments Trays



Item Code	Item Code	Sizes in mm
With Lid	Without Lid	
99-8837	DJ-4848	200x100x25MM
99-8838	DJ-4849	200x150x40MM
99-8839	DJ-4850	230x150x50MM
99-8840	DJ-4851	250x140x50MM
99-8841	DJ-4852	250x200x50MM
99-8842	DJ-4853	300x200x50MM
99-8843	DJ-4854	350x250x50MM
99-8844	DJ-4855	350x300x50MM
99-8845	DJ-4856	400x300x50MM
99-8846	DJ-4857	450x300x50MM

Instruments Box



Item Code	Sizes in mm
99-8858	150x80x40MM
99-8859	170x70x30MM
99-8860	180x80x40MM
99-8861	200x100x40MM
99-8862	230x120x50MM
99-8863	250x125x50MM
99-8864	300x150x50MM
99-8865	300x200x50MM
99-8866	350x175x65MM
99-8867	350x180x60MM

Kidney Tray (With Lid)

Item Code	Sizes in mm
99-8873	150x30 6ml
99-8874	180x33 7ml
99-8875	200x35 8ml
99-8876	250x40 10ml
99-8877	280x45 11.5ml



Kidney Tray (Without Lid)

Item Code	Sizes in mm
99-8878	150x30 6ml
99-8879	180x33 7ml
99-8880	200x35 8ml
99-8881	250x40 10ml
99-8882	280x45 11.5ml



Kidney Tray (Shallow Type)

Item Code	Sizes in mm
99-8883	170
99-8884	220
99-8885	240



Kidney Tray (With Lid)

Item Code	Sizes in mm
99-8868	250x75x40MM
99-8869	230x80x50MM
99-8870	300x80x40MM
99-8871	280x90x45MM
99-8872	440x100x50MM

Needle Cases



Item Code	Sizes in mm
99-8886	150x90x15
99-886A	210x125x17

Item Code	Sizes in mm
99-8889	50x30x8
99-8890	65x45x8
99-8891	80x30x8

Item Code	Sizes in mm
99-886C	65x15

Test Tube Stand



Item Code	Sizes in mm
99-8887	65x15


Endodontic Box




Item Code	Sizes in mm
99-8888	65x15

Dental Box

Item Code	Sizes in mm
99-8897	250x160x40




Scaler Trays



Item Code	Sizes in mm
99-8892	200x100
99-8893	205x115
99-8894	210x150
99-8895	220x130
99-8896	225x100

Bowls



- (A) Wash Basin
(B) Lotion Bowls
(C) Bowl
(D) Gallipots
(E) Iodine Cup
(Available in all sizes)

Dressing Drum




Item Code	Sizes in mm
99-8898	80x80
99-8899	120x120
99-8900	130x290
99-8901	150x140
99-8902	160x120
99-8903	180x180
99-8904	200x190
99-8905	240x240
99-8906	290x240
99-8907	340x240
99-8908	390x200
99-8909	390x290
99-8900	390x390

Bed Pan



Item Code	Perfection
99-8911	(Adult) Perfection
99-8912	(Child) Perfection



Item Code	Slipper
99-8913	(Adult) Slipper
99-8914	(Child) Slipper

Urinals

Item Code
99-8835

Male Urinal Standard Size



Item Code
99-8837

Urinal Jug Standard Size

Item Code
99-8836

Femal Urinal Standard Size





Irrigators (Deuche Can)

Item Code	Sizes in mm
99-8917	1.0
99-8918	1.5
99-8919	2.0



Irrigators (Deuche Can)

Item Code	Sizes in mm
99-8906	250
99-8907	500
99-8908	1000

Holding Pins



Item Code

99-8920

Holding Pins For Instruments



Item Code

99-820A

Holding Pins For Partition Sheets

Stainless Steel Tag



Item Code	Sizes in mm
99-8923	50x25mm
99-8924	75x30 mm

Item Code	Sizes in mm
DJ-4932	90x50 mm
DJ-4933	90x35 mm

Spirit Lamp




Item Code	Sizes in mm
99-8921	Single
99-8922	Double

Silicone Pin Mat



Item Code	Sizes in mm
99-8928	500x300
99-8929	480x250



Universal Silicon Holders



Item Code	Sizes in mm
99-8950	250x35x10



Item Code	Sizes in mm
99-8951	250x35x10

This image shows a full page of blank graph paper. The grid consists of small, uniform squares formed by thin gray lines. There are no margins, text, or other markings on the page.



Digi-med

Inhaber: Patrick Mohr

Kreutzerstrasse 1

78573 Wurmlingen

Telefon: +49-(0)7461-9101172

Fax: +49-(0)7461-9101172

info@digimed-de

www.digimed-de

Rev.: 2024